

Product datasheet for **SC211897**

TTC22 (NM_017904) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TTC22 (NM_017904) Human 3' UTR Clone
Symbol:	TTC22
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017904
Insert Size:	1025 bp
Insert Sequence:	>SC211897 3'UTR clone of NM_017904 The sequence shown below is from the reference sequence of NM_017904. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
CCTCCCAAAGTGCTGGAATTCAGGTGTGAGCCACTACACCCAGCCCTGTAAGGTTTTTGAGTGAAAGA
ATGCAGGAGCCTCCACTCACTGCAGGAATCCCTTCTACAATTCCTGGTCAGTGAGCCCCAGCCTCTTC
TGGATTGGCAGGGTGGAATCAGACTGACCAGAGTTGAAATCTGGCTTCTCCGCCTGTGGAACCTTGCA
GCAAGCCACTTACCCCCACAACTTCAGTTCCCAATCTGTAAGATATAAATAACAACACTTTCCTTAA
AATGTTGTCAACAGACAGTGAAGTAGCATCTAAAGCACAGTGCCTTGGACATAGGGGGTGCTTTTCT
GAAAGCCAGTAGAGGTCAATTAATATAGTGGATCAGAGCTCTGGCTCTGAACACAGACAGCCTTGTTTC
AAATCCCAACACCACGTGCTACGGGACCTTTGGCAAAATACTCTTGCCAGATTCTGAGCCTCAATT
TCTTTATCTGTAGACTGGAGTTGTTATGAGAATTACATTAGATTTGTATGTAAACTCAGCATCAAGC
ACACAGTAAGATTTATAAGAGACTGACTGCTACTAATTTAATATTAGGCATATTTTGAATTATCTCA
CTGTCTCTTGATGCAGTCCACTGCATAGTCATTAGCTCTGCAAGTTACACACTTAGGTGAGTTCTGCC
TCCTGCTAGCTTCCGCACTGGTCCAGACCTGCCTCCTAGAGCCACACAGAGCCCATGTGCTTGAGG
AAGCCGACTCGGAGCTCCAGTTTACTTCACTGCAAGCCGAAGAGTCCTCCCATCCCCTCATCTTCTCC
ACTACTATCCACATGACCTAGTGGCAAATCTACTCCAGTCTTGGTCTCAGTGGCTTAATTACTTTGCC
AGCAGGTACGAGGCAGGAGTGAAATTCCTTCACTGTCGTGGATGTGGTTATTATTCTAGATAAATAATGA
CATTAAACCATAATAGTAAATAATAAAATATCTCTATAAAATTCGCATTTTCATCAAGCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_017904.4</u>
Summary:	This gene encodes a protein with seven tetratricopeptide (TPR) repeats. Tetratricopeptide repeat containing motifs are found in a variety of proteins and may mediate protein-protein interactions and chaperone activity. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jun 2011]
Locus ID:	55001
MW:	38.9