

Product datasheet for SC211915

TMEM151A (NM_153266) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM151A (NM_153266) Human 3' UTR Clone
Symbol:	TMEM151A
Synonyms:	TMEM151
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_153266
Insert Size:	1041 bp
Insert Sequence:	>SC211915 3'UTR clone of NM_153266 The sequence shown below is from the reference sequence of NM_153266. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TGCCAGGGGGATGGCAGGGTGTCTCTGACACCCCCACGGCCCCCAGAGTGGCCCCCTCCCCACCAT
TCCACCATGGGCTTAGATGCCCCAGTGATTGTTGTCCAAACAGGCGGGAAAAACAGACCAGCCACACAC
AAGGGGCAGGGGTGAGGGTGGGGTGGGGTCTTAAACAGACTAAAATGCAGTTACCTGTGGTCATTT
TGGAGGAAGTGAAGCCTAAGAATCGCCCCAGCATGGCTTCATCCCCCAAGATGGCCGATGATCTGGC
CTGAAGGCCTCTCACAAGGGCAGGAGAAGAAGCTTCCAGAAAGGCTGGGATGGAGGAGGAGTCCTGCG
GCCAGCTGTTGCCTGCAGCAGAGGCTGCCTGTGGACACTTCCCTACCGGGGCTGAGGCCATGCTGGGAG
CCCCCCCAGTGACGCAGAAGCACAAATTGGAGGCTTGGGGCCATGCCAGCCAGGTTGGCTTCCCTGGG
GCTGAAGAGCTGGCTGTGGTGGCACCATGGCCTTGGCCCTGCTACTACCCCAACATTCCCATCCTCCCC
TGATAGGCCCCCTCCTCCTTGCCTGCCCTCTCTCCCGTGCTCAGGGCCAGTCACTCCCCCATGACAG
CTACCCCTACAACCTGGGGTTCCTTCTCAGGGCCTCACCTGGTCCCCACATCCCCTGGCCTCTGCTCCT
GCCTCCCCAAGTCCATGTCTGCCTCCATCAGCCACACACAGCAGGGACGCCCTTCTTGGCCTCCCG
GACTTATCTACTTCTCACTTTCTGATGATCCATGGGGCCAGGCCTAGAGGGCATCCGAGGGGCTGGG
CCTTGGCCACAGCCTCTTGGATGATGTCTTGTCCACCACGCAGGCCCTTGCTTCTCCCATCAACGAG
TCTTGGCCAGGGGACACTGACCCTCCCCCAGGGAGGAGAGGAGTTCCTTCCCTTTCAACGAGGGCT
TCCACGCCCTCTAAGGGGTGGCGGGGAAGGGAGGAGGCTGAGAAGATGCAGATTAATAAAGGTATGGA
ATGGAA
ACGCGT AAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_153266.4</u>
Locus ID:	256472
MW:	38.3