

Product datasheet for SC211920

Isochorismatase domain containing protein 1 (ISOC1) (NM_016048) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Isochorismatase domain containing protein 1 (ISOC1) (NM_016048) Human 3' UTR Clone
Symbol:	Isochorismatase domain containing protein 1
Synonyms:	CGI-111
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_016048
Insert Size:	1057 bp
Insert Sequence:	>SC211920 3'UTR clone of NM_016048

The sequence shown below is from the reference sequence of NM_016048. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CCAGAGTCGGGTCTGCTTTCCAAAGTAAGGACATTTGAAGAACTGGTATGCTACTCACTGGTGAAGGA
CAGTCAGGTGAAGGACTGTAAGCCACACAAGCTCTTCTTATCTCTACTAGAATTTAAATGTTAAGTCA
AAAACGGCTCCTTTTTGCGCCTCCTAGTGAACTTAACCAGCTAGACCATTTGAGTACCAGCATTTAG
TTACAAACGTCAAAGGCTTCCGGTGCTGCTTACCTTCCTTTTTGTTAATGTGCTTTTATTATTTAAAA
AAAATTACAATGAAGATGCCTGTTTGTCTCTACTGTGTACTCTGATCGTATCTTTCCAAAGTGCAGAC
TCTTGTGAAGTTTTCTTAAATTGTTCACTTTAAAGAAAATGACGTACCAACAATGATTTGGCTTTTATA
TTACTGTAAGATGTTATAATGTTAATGTGGATGTAGTGCTTTTACTTTACAGATTGATTGGAATAAGAT
TATTGCATATGAATTTACCCACAGGACTCTGAATCATGTTACCCACTCCCTCACAATGTTGTCCACTT
AGTGAGTTGCATTGATCTATCCGTACCAAATGATGTTGAATAATTACATATCTTTCTTGACTATACTGA
TTTCTTATTTTGGTCACTATTACTAAATCTCTGTTAATATTCTCTCTTTTAACTGAAAAGGGATGGGAT
AGAAGGGTTTGCAATGCCATATTATTGGTGGAGGGCTGTTTTAACATCTTTGAAGTATGGCTTGCTGAA
TATCTTTACCAACATCTTGAATATATATTCTAGTGTCACAAGATTTAGCAAAAAGATAAAGCTTGGGT
GGAATATCATTTTAAAAATGTTTCATGTTCTGTTCTATATTTTCTTACCTACTCTCCAAATATTGTAAATG
CAAAAAGTCTCAGTAATGATTTGGTAGTATTAATTTTGGTTCATTGTTTCTCTTCGATAAATTTATTT
TCATTAATACTTGTTAGAGGGTTTTGAAATGTTTTCAATATGTGAAATGTGAAACTGCTGTCTTTT
ATATTAAGTAATTAAGAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_016048.2</u>
Locus ID:	51015
MW:	40.5