

Product datasheet for SC211921

PPP1R11 (NM_021959) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PPP1R11 (NM_021959) Human 3' UTR Clone
Symbol:	PPP1R11
Synonyms:	CFAP255; HCG-V; HCGV; IPP3; TCTE5; TCTEX5
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_021959
Insert Size:	1057 bp
Insert Sequence:	>SC211921 3'UTR clone of NM_021959 The sequence shown below is from the reference sequence of NM_021959. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CAGCCCCCTCCAGGGCAATGCAGCACTAATCCCTCTCTCCTCCAGCATTCTGTGTCTGTCTGGCCC
TAAATGTATCCATGTGGTACTTCTCCAGCCCCCTCCTCCCTCTCTTCTGCCTGATAGAGGGAAGAGG
AAGAGGAGGACGAACAGAGATCCTGAAATTCTGACTTGCTGCTATTCCAGAACCAGCCTCCTGGGTTT
CCCCAGTCCTCATTTTTCTCCCAATACCCACCCTTCTCTCTCGAGGGATCTAGGCACCTTGGTCCCAG
TGTCTTCTTTTGTCTCACTGCCAAACTGCCTGTCTGGGATCTAGTTATCTTGCCCTGCACTCTCA
ACATGAGTAGCGAACACTTAAATTGGGTTTTCAACAGTCCCAGCTTCACTGCCAGGGTCCCAGTCAGA
TTCCAGGAATTTGCGCCCTAACTTTGCTTGTAACTCCTGGTTTAGAGCTATCCCACTAAATATTTAAT
CCTAATCTTAGTCCTTGCTGTGAGATATGAGGTCTTACAGGAGACCTCAGAGCTCCAGCCCTTCTC
CTCCTGCTAACCTTCTCACACCCTCAAGAGGAGTTAGAAAAGAGGTCTTGTCTTCTCACCTCTTAT
GGAAAATGGAATAAGAAATAATCATATCCTTTCTTCCACCCTTCTCCTGTTATTTAGGATTTCTGACA
AAGCTGGCTTGAGATTGGTCACTTAGAGCCGACTGTCTCCTCTGCCTTTTGTCTTTCAGCTTCAGAGAC
AGATCCAATATAGTCCCAGGGACCTGGGTCTCTGGGAGAGGAAGGAAGAGGGAGGCAAGAGATTG
GGGTATGTCCCCTGTAGTACACTCTTACCTCTTACTTCTAGACTTTGATTTCTCCGGCAGCCAGATG
TTCAGTTCTCTTGGCCCTCTCTACCCCTTACTGGGATCCGGTTTTTCATTTTCCGGTCTTTTGCCATA
CACAGTTACAGAGATCAGTCAAATCCATACCACCACTGAGATCTCATTTATTGCCACAGATGCACAAAA
TAAATAACCCAAATCACAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_021959.3
Summary:	This gene encodes a specific inhibitor of protein phosphatase-1 (PP1) with a differential sensitivity toward the metal-independent and metal-dependent forms of PP1. The gene is located within the major histocompatibility complex class I region on chromosome 6. [provided by RefSeq, Jul 2008]
Locus ID:	6992
MW:	39.4