

Product datasheet for SC211924

TMEM9B (NM_020644) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM9B (NM_020644) Human 3' UTR Clone
Symbol:	TMEM9B
Synonyms:	C11orf15
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020644
Insert Size:	1148 bp
Insert Sequence:	>SC211924 3'UTR clone of NM_020644

The sequence shown below is from the reference sequence of NM_020644. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTCTTTGACCGGCATGTTGTCCTCAGCTAATTGGGAATTGAATTCAAGGTGACTAGAAAGAAACAGGCA
GACAACTGGAAAGAACTGACTGGGTTTTGCTGGGTTTCATTTAATACCTTGTTGATTTACCAACTGT
TGCTGGAAGATTCAAACTGGAAGCAAACTTGCTTGATTTTTTTCTTGTTAACGTAATAATAGAG
ACATTTTTAAAGCACACAGCTCAAAGTCAGCCAATAAGTCTTTCTATTTGTGACTTTTACTAATAA
AAATAAATCTGCCTGTAAATTATCTTGAAGTCCTTACCTGGAACAAGCACTCTCTTTTACCACATA
GTTTTAATTGACTTTCAAGATAATTTTCAGGGTTTTGTTGTTGTTGTTGTTGTTGTTGTTGTTG
TGGGAGAGGGGAGGGATGCCTGGGAAGTGGTTAACAACCTTTTTCAAGTCACTTTACTAAACAACTTT
TGTAATAGACCTTACCTTCTATTTTCGAGTTTCATTTATATTTTGAGTGATAGCCAGCCTCATCAAG
AGCTGACTTACTCATTTGACTTTTGCACTGACTGTATTATCTGGGTATCTGCTGTGCTGCACTTCATG
GTAAACGGGATCTAAATGCCTGGTGGCTTTTCAAAAAAGCAGATTTTCTCATGTACTGTGATGTCT
GATGCAATGCATCCTAGAACAACTGGCCATTTGCTAGTTTACTCTAAAGACTAAACATAGTCTTGGTG
TGTGTGGTCTTACTCATCTTCTAGTACCTTTAAGGACAAATCCTAAGGACTTGGACACTTGAATAAAG
AAATTTTATTTTAAACCAAGCCTCCCTGGATTGATAATATATACACATTTGTCAGCATTTCCGGTCGT
GGTGAGAGGCAGCTGTTTGAGCTCCAATGTGTGCACTTTGAACTAGGGCTGGGGTTGTGGGTGCCTCT
TCTGAAAGGTCTAACCATTATTGGATAACTGGCTTTTTCTTCTATGCTCTTTTGAATGTAACAAT
AAAAATAATTTTGAACATCCATCAGTGATCTATCTATGCTCCTAGTTTTTCTCCTCCTCCTCTTT
TGCTGTATAATGAGATTGAAGATATAAAGACATTTGTACCTG
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_020644.3</u>
Summary:	Enhances production of proinflammatory cytokines induced by TNF, IL1B, and TLR ligands. Has a role in TNF activation of both the NF-kappaB and MAPK pathways.[UniProtKB/Swiss-Prot Function]
Locus ID:	56674
MW:	44.8