

## Product datasheet for SC211928

### UNC5A (NM\_133369) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** UNC5A (NM\_133369) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** UNC5A  
**Synonyms:** UNC5H1  
**ACCN:** NM\_133369  
**Insert Size:** 1042 bp  
**Insert Sequence:** >SC211928 3'UTR clone of NM\_133369  
 The sequence shown below is from the reference sequence of NM\_133369. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCTTCACAGTGTCCGAGGCTGAGTGTCTAGGCCCGCCAGGCCGACACCTACACTCTACCAGCTTTG
GCACCCACCAAGGACAGGCAGAAGCCGGACAGGGGCCCTCCCCACACCGGGGAGAGCTGCTCGGACAG
GCCCCCTCCCGGCCGAAGCTGTCCCTTAATGCTGGTCCTTCAGACCCTGCCCGAACTCCACCTCTCCA
TGGCCTGCCTAGCCAGGCTGGCACTGCCACTCACACTCGGCCCCAGGGCCAGGAGGGACAGTGCCTGG
AGCCTGGGCCAGGCCAGCCATCTGTGTGTGTATGTGCGTGTGATGCTACCTCTCTCCCGTCCCT
CTCCAGGGGCCCGCATAACACACGGCCATGCACGCACACACTGGGCCTGGGCCAGGGCCCGAGACTCC
TGCCTGAGCTGGACCTTATGCAAACATTTCTGTGCCTGCTGGGTAGGGGCACGTCTGAGGGGCCCTGCT
CCAAGCCTGCAGGACCGAGGGCCACAGCCGGACAGGGGTAGCCCCTGGATTAGGCACACGACCACCA
CACGAGCACGTGCCACGCATGCCTCGTGTGCTCATCTCACACACACCCCTCCCGGGTCACGCAGACA
CCCCCAACCACACACATCTCATGCTGTACACCTGAGGCTGCTCACGTCTCACGCCAGTGTGGTGA
CATTTGCCTCTCATATGCTGCCCTCTCCACCCACCCAGGACACCCACGGCTCTCCCTGCCCTGCC
CCTCCCCAGCCTTGAGGTGCCCTGCCCGGGCCGTGAATATGCAATGGGAGTCCCAGGCTGTAC
AGTGGTGAAGTGTGTGTGGCGTGGCGTCCCCAGGGCTGGTGGTCCCCACGCGGGCCCTGT
CATGTGAAGCTCGTGTCTGACTTTGTCTTAAGTGCAATCACGCCTTACTCTTGGCCTTATGTACACA
GCCTTGCCCGGCCCGGGGCACATAGGGGTTTTATCGGGCGTGAATGTAATAAATTATATATATA
TTGCTAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

|                        |  |
|------------------------|--|
| <b>OTI Disclaimer:</b> | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| <b>Components:</b>     | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.   |
| <b>RefSeq:</b>         | <u><a href="#">NM_133369.3</a></u>   |
| <b>Summary:</b>        | UNC5A belongs to a family of netrin-1 (MIM 601614) receptors thought to mediate the chemorepulsive effect of netrin-1 on specific axons. For more information on UNC5 proteins, see UNC5C (MIM 603610).[supplied by OMIM, Apr 2004]  |
| <b>Locus ID:</b>       | 90249  |
| <b>MW:</b>             | 37   |