

Product datasheet for SC211956

CCRK (CDK20) (NM_012119) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCRK (CDK20) (NM_012119) Human 3' UTR Clone
Symbol:	CCRK
Synonyms:	CCRK; CDCH; P42; PNQALRE
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_012119
Insert Size:	1051 bp
Insert Sequence:	>SC211956 3'UTR clone of NM_012119 The sequence shown below is from the reference sequence of NM_012119. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGATTCGGCCCTTCATCCTGGAGGGGTGAGAAGTTGGCCCTGGTCCCGTCTGCCTGCTCCTCAGGACC
ACTCAGTCCACCTGTTCCCTCGCCACCTGCCTGGCTTACCCCTCAAGGCCTCCCATGGCCACAGTGG
GCCACACCACACCTGCCCTTAGCCCTTGCAGGGTGGTCTCGAGGCAGAGGTCATGTTCCAGCC
AAGAGTATGAGAACATCCAGTCGAGCAGAGGAGATTTCATGGCCTGTGCTCGGTGAGCCTTACCTTCTGT
GTGCTACTGACGTACCCATCAGGACAGTGAAGTCTGCTGCCAGTCAAGGCCTGCATATGCAGAATGACG
ATGCTGCTTGGTGTGCTCTCCCGAGTGTGCTCCTGGTCAAGGAGAAGTGCAGAGAGTAAGGTGT
CCTTATGTTGAAACTCAAGTGGAAAGGAAGATTTGGTTTGGTTTTATTCTCAGAGCCATTAACACTAG
TTCAGTATGTGAGATATAGATTCTAAAACTCAGGTGGCTCTGCCTTATGTCTGTTCTCCTTCATTT
CTCTCAAGGAAATGGCTAAGGTGGCATTGTCTCATGGCTCTCGTTTTGGGGTCATGGGGAGGGTAGC
ACCAGGCATAGCCACTTTGCCCTGAGGGACTCCTGTGTACTTACATCACTGAGCACTCATTTAGAAG
TGAGGGAGACAGAAGTCTAGGCCAGGGATGGCTCCAGTTGGGGATCCAGCAGGAGACCCTCTGCACAT
GAGGCTGGTTTACCAACATCTACTCCCTCAGGATGAGCGTGAGCCAGAAGCAGCTGTGATTTAAGGAA
ACAAGCGTTCCTGGAATTAATTTATAAATTTAATAAATCCAATATAATCCCAGCTAGTGCTTTTTCT
TATTATAATTTGATAAGGTGATTATAAAAGATACATGGAAGGAAGTGAACCAGATGCAGAAGAGGAAA
TGATGGAAGGACTTACGGTATCAGATACCAATATTTAAAAGTTTGTATAATAATAAAGAGTATGATTGT
GGTTCAAGGATAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_012119.5</u>
Summary:	The protein encoded by this gene contains a kinase domain most closely related to the cyclin-dependent protein kinases. The encoded kinase may activate cyclin-dependent kinase 2 and is involved in cell growth. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Dec 2009]
Locus ID:	23552
MW:	39.1