

Product datasheet for SC211985

Centaurin alpha 1 (ADAP1) (NM_006869) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Centaurin alpha 1 (ADAP1) (NM_006869) Human 3' UTR Clone
Symbol:	Centaurin alpha 1
Synonyms:	CENTA1; GCS1L; p42IP4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006869
Insert Size:	1051 bp
Insert Sequence:	>SC211985 3'UTR clone of NM_006869 The sequence shown below is from the reference sequence of NM_006869. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTGGAGGCGCACTTCAAGCATAAACCTACGAGTGGGCTGGAGGACCACGGACATTGGACTCACTGT
GGCTGGACGGAGGGGACCCGTGGATGGGGGGCTCTGGCGTCTGAGGCCACCTGGCCCCACCTGCTCC
TCAGGGCAGCCCGCGCGGCCAGGTAGGCGCCGAGCTTCAGCTTCAGGATGCTTCTCTGGAACCTCAA
GGCAGGCAGCCCAGGCCCTGGGCTGATCTCTAACCCCGTCATGCTGCTGCTGACCACACCCAGCCAAC
CTGCCCCCTCCCTGACCCGGGGCCCCCTTTCTCCAGGGCCAGTCTGGCTCCCGAGCTCAGACACAGCC
CCCAGGAGCCCGCCACCACGGCTGTGGGAACCTCCAGGCCTCAGCGCAGAGCCCTGCTTGGGGCCAGCA
CGGTGGGCGAGTGGCGCCGCTGGGGACTGGGCCCTGCTCGCATGCCGCCCGCCCTCCCCCACCTCC
ACGACTATTTATTGAGCGCCTGTTGTGTGTACGGGGCTATGAGGGCCGTGGGGTGTGGGTGGATTA
TCCACACAGGTCCCGGCCCTGCCCGGGCTGGAGTTGCCACAGCCTGTGCTCCTGGTCTCACCTGGAG
GGGCCAGCAGGCTGCCGTCCACACACGTGGCCTCTGCGCCAGCACGGTGCTCGCCGACAGTGGTGT
CTGAACCTTTGGGACGAGGGCCTGGGCCGCGGTGAGGCCACAGAGGACAGGAGTGGCCCTGGGGTCC
CGGGCACTGTGCGCTTGTGTCAGGCGGCCAGCCGTGATTTATTTTACCTATCTCCTTCTGTCA
AGGCAGGCCGGGCTCCAGGGCTCCCTTGCCTGGGGCATGTAGGTGGGGGAGGCGTCTGCAGGTACCT
GGGGGGCCAGCCCTCTCCAGCCTTGCTGGCTGAGCTGTGTTCCAGGGAGCCCTGGACAAGCCCTC
ATAGGCAGGGAGGGGTTTCCGAGGCCAGGCATCCGCCGCCCGTGTGCGATCCTGGAATAAAATGTGG
CTCTGGCATGGTTCCA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006869.4</u>
Summary:	GTPase-activating protein for the ADP ribosylation factor family (Probable). Binds phosphatidylinositol 3,4,5-trisphosphate (PtdInsP3) and inositol 1,3,4,5-tetrakisphosphate (InsP4).[UniProtKB/Swiss-Prot Function]
Locus ID:	11033
MW:	37