

Product datasheet for **SC211991**

CCDC191 (NM_020817) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC191 (NM_020817) Human 3' UTR Clone
Symbol:	CCDC191
Synonyms:	KIAA1407
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020817
Insert Size:	1048 bp
Insert Sequence:	>SC211991 3'UTR clone of NM_020817 The sequence shown below is from the reference sequence of NM_020817. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TTGAGTAAGACAAGTCTGGTGAACGAATAGCCAGATTAGTGAATATTTAGGACAGGAAGCTCGACTA
GTCTGTATGTGTATTCTCCACCCACCCACACATACATACACCTGGTTTATTCTCTACGAGC
TACTTCTTAGATCAAAGGAATTATTTGAGAATTTTATTCCTTCTTTGACTTTTCTGGATTGCTAG
TCATTCATATAGCATTTATCACAGTGGGATCATCCAAAGGAACATTCTAGAGATCGTCTGTTCAAATA
TATCCTATATATTTTAGCATAGTAACTTTAAAAAGGACAAGTGTTCCTGACTGAGCCTGTGAAGA
AAATGTGTCCGTATAGCACGACAGGAGTTGCTGATGTGACCTTAATGAAACAGAAATAAACAGGGGTGA
GTGCTCGGAAAGGAAGGGCGGGCCAAGGCACTCAGCCTTCTAGTGCAGTAGCATCAGTGACACGACCCC
CACCAGACCTCGTCCACTTTCACCAACCCTGAAGGCAGACCATGTATGAAAAGGTTCTCTACACAGCTG
CTGCAGAGGACAGCAGCTTATCCTGATTTCAGAGAACTGAGGCACCAGAGATGCTTCTCTCCGTGCA
CCTTGATAAGCATGAAGGACATTCAAAAGGTAATAGCTCTATCATGACAGTCCAGCAATTCATTAAT
TATCTGTTTCATGCTCCACCTGCAACTGCTGAATAATACACTCAGTTTGTCTTCTAGTACTCAGCTCAG
TTTGATCTAGATATTGAGATACCTATCTGAATTCAAAGAATGTTTTCATCTTTAAGATGCTATAAAGA
CTGATAAAATGCTCTCTTCTATGAATTGTGTTAGAGAATTACCTGCTAGAAAAATAATTTTTCATTATA
AAATACACATTTATTATATAACAAATTCAAGTAGTCCATAAATGCATATGGTAAAATGAGAAAGTCTG
TGCAGTCTGAGATGATTACTCTATTTGGTGTGCTGAAAATTTGAGTTGACAAAATAAAATTCTGTGG
TTTTTAAAGTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_020817.2</u>
Locus ID:	57577
MW:	40.4