

Product datasheet for SC212008

ZBTB9 (NM_152735) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZBTB9 (NM_152735) Human 3' UTR Clone
Symbol:	ZBTB9
Synonyms:	ZNF919
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152735
Insert Size:	1052 bp
Insert Sequence:	>SC212008 3'UTR clone of NM_152735 The sequence shown below is from the reference sequence of NM_152735. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TGGGCCACTGCCACTGGACTTACAAGTGACTGTGCTGAGGCTATACACTAGCTTCTAGAACAAGATAACC
ACTGCTGCTGATGGATACTTTTCCCTCACTGCCATGGCACACCAGTCATGGATCTTGTAAATCATGCCAA
GAGAATAGATACATTATGGACCTCTTGTTCTTAGATATGGGCCTCTCAGCCTGGCAGATGTGAAACTC
AAATTTCTCGTCCCACTCCAGGTTTTGGCTAGCCAACCCTGCAGGAAAGTGGTTTATAGGCCATTGATA
CTTAAGTTGATCACTTGCCCATGGTGGACATTTTGTGGTGGTGATGTCCATTAAGGAAACCAGATTTT
CAATTATTTAGTGAGAGAAGAGTTAGAGCAAAAGACAGTGGTAAATGTTTTATTCCGTCTCCATGAGGA
ATTGAAGGAGTTGGTCTCCACCTAGAGATACATTTGATTTACAGCTTAAGTAATTCAGAGGCTAAGCTC
TAAGCTTTTTTCTCTCATTGCTGGAATGATTTAAGCAGAAGTCCTTTTGTGTACTTTTAAATTTGTATC
TTTCCAGGAGCCCTCAGATTGTACCTTGCTTTCTACCAATAGACACCTTCCGACACTTTTTTAATG
TTGTAGCTGAGCACTTTAACAAGTTGAGCATTCCATGTTTCATTCTTAGAACCTTCTTTAATAGAGGGT
CTTCCCTCAACAGCCTGTGCCTCTGGTCTACCTTTGACCACCACTGATACTAATATATTGGTCACAAT
GACTGGAATGTGACTAGTGATCTCAGGAGATGGCACTGTCCTAAAGTGCTGTCAGGGTGGCACCCTGC
TCTCTGAACAACCTTACCTTGGTGAGAGGACTCAGGTTTGGGACAGCACAAAGCTGAAGGCTGGAGAGTA
ACTTGATAGTAGGACCATACCTCTTCCCTTCCCATCCACCCACATATGATAGACAGCCCTCTGTTG
AGATATGGAGGGGACAGATACTGGAATCGGGGTGGGACTTGCACTTAAATTTTTTAATAAACT
GTGCCCTGAAACCTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_152735.4
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	221504
MW:	40.5