

Product datasheet for SC212021

UNG (NM_080911) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: UNG (NM_080911) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: UNG
Synonyms: DGU; HIGM4; HIGM5; UDG; UNG1; UNG2; UNG15
ACCN: NM_080911
Insert Size: 1053 bp
Insert Sequence: >SC212021 3'UTR clone of NM_080911
 The sequence shown below is from the reference sequence of NM_080911. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGAAGCCCATTTGACTGGAAGGAGCTGTGATCATCAGCTGAGGGGTGCCTTTGAGAAGCTGCTGTAA
CGTATTTGCCAGTTACGAAGTCCACTGAAAATTTTCTATTAATTCTTAAGTACTCTGCATAAGGGGG
AAAAGCTTCCAGAAAGCAGCCATGAACCAGGCTGTCCAGGAATGGCAGCTGTATCCAACCACAAAC
AAAGGCTACCCTTTGACCAAATGCTTTTCTCTGCAACATGGCTTCGGCCTAAAATATGCAGAAGACAGA
TGAGGTCAAATACTCAGTTGGCTCTCTTTATCTCCCTTGCTTTATGGTAAAACAGGGGAGATGTGCAC
CTTTCAGGCACAGCCCTAGTTTGGCGCCTGCTGCTCCTTGGTTTTGCCTGGTTAGACTTTCAGTGACAG
ATGTTGGGGTGTTTTTGCTTAGAAAGTCCCCTTGCTCAGCCTTGCCAGGGCAGGCATGCCAGTCTCTG
CCAGTTCCACTGCCCCCTTGATCTTTGAAGGAGTCCAGGCCCTCGCAGCATAAGGATGTTTTGCAA
CTTTCAGAATCTGGCCAGAAATTAGGGCTCAATTTCTGATTGTAGTAGAGTTAAGATTGCTGTGA
GCTTTATCAGATAAGAGACCGAGAGAAGTAAGCTGGGTCTTGTTATTCCTTGGGTGTTGGTGAATAAG
CAGTGGAATTTGAACAAGGAAGAGGAGAAAAGGAATTTGTCTTTATGGGGTGGGGTATTTCTCCT
AGGGTTATGTCCAGTTGGGGTTTTTAAGGCAGCACAGACTGCCAAGTACTGTTTTTTTTAACCGACTGA
AATCACTTTGGGATATTTTTCTGCAACACTGAAAAGTTTTAGTTTTTTAAGAAGTACTCATGCAGAT
ATATATATATATTTTTCCAGTCTTTTTTTAAGAGACGGTCTTTATTTGGGTCTGCACCTCCATCCT
TGATCTTGTTAGCAATGCTGTTTTGCTGTTAGTCGGTTAGAGTTGGCTCTACGCGAGGTTTGTAAAT
AAAAGTTGTTAAAAGTT
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_080911.3</u>
Summary:	This gene encodes one of several uracil-DNA glycosylases. One important function of uracil-DNA glycosylases is to prevent mutagenesis by eliminating uracil from DNA molecules by cleaving the N-glycosylic bond and initiating the base-excision repair (BER) pathway. Uracil bases occur from cytosine deamination or misincorporation of dUMP residues. Alternative promoter usage and splicing of this gene leads to two different isoforms: the mitochondrial UNG1 and the nuclear UNG2. The UNG2 term was used as a previous symbol for the CCNO gene (GeneID 10309), which has been confused with this gene, in the literature and some databases. [provided by RefSeq, Nov 2010]
Locus ID:	7374
MW:	40.3