

Product datasheet for **SC212037**

TMEM71 (NM_144649) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM71 (NM_144649) Human 3' UTR Clone
Symbol:	TMEM71
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144649
Insert Size:	1056 bp
Insert Sequence:	>SC212037 3'UTR clone of NM_144649 The sequence shown below is from the reference sequence of NM_144649. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACTGCCTGTGCTCGGTTTGTCAAAATTTGCAACCATTTAGGAATGCCTTCGATGAATGTCCTCCATCT
GAATATCTGGAATTGTCCAACCTGCAGTCTACTTGAATCAAGTGTTTTATTGGAAGGGAGTAAGCGA
GTAATGGAGAAAAAGCCATTTTAGTTTGACTATGTGATTTTAAATGATCTCAGTTTTCCATCAAAAT
TATAATATGCTCATGAAAATAATATTAATTTGCCTTCCCTTTGCAACACCGGCAGTTGAAAGGAAAAG
GACGGGGAATGTGATGGAAAAGAGACCGCCTGGAATAAATGTCCCCTATGATTCTTTAAGGCAGTGTT
TCTCGAGCTTGAATTTTCATTAGGAAATCTGTGAGGAGCTTGTAAACAGATTTCTGGGTCTGCCACAT
GCACCTATCTCTTGCTGAATTGCTTTAATAGAATAATGAGAGCAAGTTTGTCTAACTAATACCAACCTG
ACAACCTGAATAACAATAAATGCAATTTGTACATAAAATATAATGCTGCAAAAGTTTGTCAATTCACCTC
AGTGGAGTGACTTGATATTAGGTGGTAACCGTAGATGATGGTTAATGAAAATGGACAGGAAAGAAGC
ATTTTCTGAAAGTTATATTCTTTTGAACACGTTCTAAACCAAGTTTTTAATCTTCTTGGGGCTCGTAA
TTACCTTTCACTTTAATGTCACCTAAAGATATAACACAGAAAAATGCCTTGAGGGCAAAATATAGGCAA
AACACCAATGCGCTTTCAAATGCATGAAAATGGTGCAGTTGTACCCTTGAGCCTTGACTCAAGGGCTGT
AGATGTTCCCTTTCCACCCCCACACTTGGTGCCTGTTACAAAGCAAATATGGCCTGTAATTCAAATT
TGTTCATGTGATACTCTCTGAGTAAAACTCATACATGCAGAAAATTGTCTTTGCTCGAAATGATAAT
GCCAAATATACTTTATATATAATTTGCATTTAGTACATTTTGGTTAAAAAATAAACTAATAAATAA
GTGAAGTCATCAGCAATTTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_144649.3</u>
Locus ID:	137835
MW:	40.2