

Product datasheet for SC212044

C8orf44 (NM_019607) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C8orf44 (NM_019607) Human 3' UTR Clone
Symbol:	C8orf44
Synonyms:	FLJ11267
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_019607
Insert Size:	1234 bp
Insert Sequence:	>SC212044 3'UTR clone of NM_019607 The sequence shown below is from the reference sequence of NM_019607. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CTTCTGTTCTCTAACCTTGAAGGTTTGAACAATATTTCTCAGAGGATAGCCTTTCACTTATTCATCT
GTCCAGCATGACTCATCCCCGGGAGTGTTGAGTAAGTGAAATTTTGCTGTATTCATGTTTTGTGACTT
ATAAAATAGGATGATAAGGAGAGAACATGAACCTCTGGAGTCAGACCTGTTACCTCGGACATGATACTCT
TAGCTTTGTCAATTTAGTATTTGAGTAATTTGGGCAAGCTAACATCTCTGGTCGTTCTCATCTGTAAAA
TGAGAATAAATGAAACCCACTAACCCAGAAATGGTATGAAAATGAAATGTGGCAGAAAAAAATGAAAGT
GAATAGTATCACCACTGACACACAAGCACTAAAGGCCCTTCTGTCTCCATCAGGTATGGATTTGGGGC
AACATTTGGCCAGGTCTTGTATCTTTCTGTTTCTATCTATTCTGTCTAATTCAGTGCTTTGTTTACAAC
GAATGTCTTACAAATGCTGACTGAACACTAGCATACCTGCATGAACAACAGGTAAATAAATTTAGATG
TGTTTTAATGTTTATTAATCTATCCTGTCAGAGAAGAACTGCCAGTTATAGATAAATATGATGCCAGGT
CAGGGCTGAAGAGTTGGGCAGGTTGTTATCTGCATGGGGTCACTAGGTTCCAGTGAGAGGTGGGGGCT
AAGCTCTCACCCGCTCTGCAGCCACCTGGCACCCGGTTTCAGTTTCCTGAAAGGGAGCCTTCTACTTGC
TGACGACTGCCTCATCTCTCTGAGGTTTCCTCTGATAAACAATTTCTCTGCTTTTTTTTTTTAATTA
TGAAATACTTTAAATGTAAAGGAGATAATGTGACACACATTTACCCATAATTGAGATTGCCATAATTGC
TATAACTTTTTCAAATTTTGACTTAATTTGAGACTTTTAGAAAAATTTGGCCAGGTGTGATGGCTCATG
CCTGTAATCCCAGCACTTTGGGAGGCCAAGGTGGGTAATTAATTTGAACCCAGCAGTTTCGAGACTTGCC
TGGGCAACATAGTGAGACCTTGCTCTACTGAAAACAAACTAAAAAAAAAATTAGCCAACCATGGTGGT
GCATGCCTGTAGTCTAGCTACTCGGAAGCTGAGGCAAGAGGATAGCCTGAGCCAGGGGAGTTGGGG
CTGAAGCCTGGGTGACAGAGGAAGACCTTGCTTCAAAAAACAAAAAAAAAAAAAAAAAAAA
ACGCGTAAGCGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_019607.2</u>
Locus ID:	56260
MW:	47.3