

## Product datasheet for **SC212090**

### RASL11B (NM\_023940) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** RASL11B (NM\_023940) Human 3' UTR Clone  
**Symbol:** RASL11B  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pMirTarget (PS100062)  
**ACCN:** NM\_023940  
**Insert Size:** 1060 bp

**Insert Sequence:** >SC212090 3'UTR clone of NM\_023940  
 The sequence shown below is from the reference sequence of NM\_023940. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
GCCAAAGTGAGGACTGTCACCTCCGTCTTGAAGCAGGAGGAGCACTCAAGGGGGTTTGGTCTTCCCAGGA
AGAGGGCCTGAGGTTCTCCTAGTGCAGGAACGTTGAATATTGGCAATGATTCCTGGTTCCAGAAAGGGC
TGGAGCAGAAGGGCCAAGAGGGCCTGTGGAAGTCTACAGAAAAGGAAGTGTGTCTGAGCAGGGGGA
CAGGATTGATGAGGCTTGAAAGAGCCCACTGAGCCACTCTCTGAATATGTGAAATGACTCTGTGTCTT
TTCCTTTAGAGTGGGGAGGGGCATAATCGTTTCGGTTTCTGCATTCAACTACCTTGAAATGGTGGTC
CGTTGCAGTTTACCAAATGTATTGATGTGATTTACAGTGGGAATGAAAGAACAGATTAAGCATTGATA
GGCTTTCTTTTTCCCCCTTTTTCTTTTTCTTTTTCTTTTTCTTTTTCTTTTTCTTTTTCCCAAGAGG
AAGAATTGCTTTTCTTTACCAGTATGATTTGTAAGTCTGGAACACTATTCTTACAGAATGCCTTTCTA
ACCTGAAGGATACCCAGATTTCTTTCTTTATGTACAAGATGGAAAATCCCTACCCCTCAAAAACAGGT
TGAGTTTATGGGCCAGAATATTTGTATACCAGACATTGGTAAGCTCTCATGGTTTACAGGAAGGCTCT
GGTTCCCCCTCACAGGATGAGCATGCTCAGTTGGGCGTGTGTGGAGGAGCTGTGAACTTCACGCTCTTC
AACCCAGACTTTGACTGTATAGTTTTTTGTTTTTTAACTGTTCCATCAGGAAAAAATGTGTCAGTT
GATTTTGGTGTGACTTGTGTAATGGTTTATGGCAATGAGCTGGGTTGCTTCTAGGCCTGGCTGAGT
TGTGCCTAAGGGTGGCTGAAATACTAAACACTTATCTTACAGCAAGTGAACAGGGGCTACCTGCCACA
TCCCTCCACAGATGCACCTTTAAAAAGCCACTCATGCTTTGGCTTAACTGTAATTAATTTATTTATG
TACAATAAATCGCATTGAAAAAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_023940.3</a></u>
<b>Summary:</b>	RASL11B is a member of the small GTPase protein family with a high degree of similarity to RAS (see HRAS, MIM 190020) proteins.[supplied by OMIM, Nov 2008]
<b>Locus ID:</b>	65997
<b>MW:</b>	40