

Product datasheet for **SC212155**

NPTN (NM_001161364) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: NPTN (NM_001161364) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: NPTN
Synonyms: GP55; GP65; np55; np65; SDFR1; SDR1
ACCN: NM_001161364
Insert Size: 1079 bp
Insert Sequence: >SC212155 3'UTR clone of NM_001161364
The sequence shown below is from the reference sequence of NM_001161364. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAAACTTGCAGAGAAACACAAATTAAGTACTGCTTACAATATCTTTAGGTTCTGAAACTGGTGG
CAACATGACCTGCTAAAATTTCTGCTTGGACCTCTTTGGTTCTCTCCCCTTCAAGTGAGCAACCA
CAATGACTGTCTAAAGCATGCCTTATTTAGCCTCTCCTGTAAGGGTGATCTAGCCAGGTACATTTTAAA
CAATGCTTCAGTGTAGAAGGTGTAAGTATTTTGGGCTTGATGTGCTGTGAATGTTGCTTTTTTTTTTC
CTTTGTTAAAATATTTAAATAGAAGTGAAAAGTCTCTGAGGATCAGATCATGCATGCGCCATTTTTT
ACTTAATGCAGCTGTTAAATTGGCAAAGCTCTAAAATGCACTGCTGCCATCTAGTGATACACTTTTGTA
AAGTACAGCAAAACCTACAGATATACAGTATATAAATATATATATATATATATTTTATTTTTGGGG
GTGGGAGAAAATCCAAAATAAAGTAAATGCTTGTTCATTTTTAAGCTGCTGATATTCATTCCTTATGT
ATGTTGTCAGATGAGGAAATGTGCAGTCTGGTACATAAAGATGAGTAATATAAACTGAAATCTATAA
TTTTAAGGGCTTAACCTGTGACTTTAATAAGCTGGAACAGTCCACTGAATGGGTATAATGAATTGCAGT
ATATACGTATGATTGCTTTTTAAGTGATTATCTTTCTTCTGTTAAGTCATGTAATTCATAAATCCTT
TTGCACTGATGTGTTGAACCTTATTCTGTACATTCATTCAATCAAGGCAAACCTTTTATAATTTTTCTT
TTGTTTCCAATGACCTTGAAATGTTATAGCATGGTAATATTCTATGCAACTATAGTTATACTTTTTGGT
TTGACACTGTATTTTTTACATTGATTTACTGGTTGATGATAGATTTTATAACCTAACGGTTCTCATGC
GGTGCCTAATTGTAGATGCATGTACTTGTGTGTTTGTGAATTTGAAGTGCAATGATGTATAAAAA
AGTGGATTACCTGTTTTTAAAAATAAAACATTGATAAAAGGTG
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001161364.2</u>
Summary:	This gene encodes a type I transmembrane protein belonging to the Ig superfamily. The protein is believed to be involved in cell-cell interactions or cell-substrate interactions. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2009]
Locus ID:	27020
MW:	41.8