

Product datasheet for SC212161

NPTN (NM_012428) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: NPTN (NM_012428) Human 3' UTR Clone
Symbol: NPTN
Synonyms: GP55; GP65; np55; np65; SDFR1; SDR1
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_012428
Insert Size: 1079 bp
Insert Sequence: >SC212161 3'UTR clone of NM_012428

The sequence shown below is from the reference sequence of NM_012428. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAAACTTGCAGAGAAACACAAATTAAGTACTGCTTACAATATCTTTAGTTCTGAAACTGGTGG
CAACATGACCTGCTAAAATTTCTGCTTGGACCTCTTTGGTTCTCTCCCTTTCAAGTGAGCAACACCA
CAATGACTGTCTAAAGCATGCCTTATTTAGCCTCTCCTGTAAGGGTGATCTAGCCAGGTACATTTTAAA
CAATGCTTCAGTGTAGAAGGTGTAAGTAACTATTTGGGCTTGATGTGCTGTGAATGTTGCTTTTTTTTTTC
CTTTGTTAAAATATTTAAATAGAAAGTGAAAAGGTCCTCTGAGGATCAGATCATGCATGCGCCATTTTTT
ACTTAATGCAGCTGTTAAATTGGCAAAGCTCTAAAATGCACTGCTGCCATCTAGTGATACACTTTTGTA
AAGTACAGCAAAACCTACAGATATATACAGTATATAAATATATATATATATATATTTTTATTTTTGGGG
GTGGGAGAAAATCCAAAATAAAGTAAATGCTTGTTCATTTTTAAGCTGCTGATATTCATTCTTATTGT
ATGTTGTCAGATGAGGAAATGTGCAGTCTGGTACATAAAGATGAGTAATAAACTGAAATCTATAA
TTTTAAGGGCTTAACCTGTGACTTTAATAAGCTGGAACAGTCCACTGAATGGGTATAATGAATTGCAGT
ATATACGTATGATTGCTTTTTAAGTGATTATCTTTCTCTGTTAAGTCATGTAATTCATAAATCCTT
TTGCACTGATGTGTTGAACCTTATCTTGTACATTCATTCAATCAAGGCAAACTTTATAATTTTTCTT
TTGTTTCCAATGACCTTGAAATGTTATAGCATGGTAATATTCTATGCAACTATAGTTATACTTTTTGGT
TTGACACTGTATTTTTTACATTGATTTACTGGTTGATGATAGATTTTATAACCTAACGGTTCTCATGC
GGTGCGTAATTGTAGATGCATGTAAGTGTGTTTGTGTAATTATTGAAGTGCAATGATGTATAAAAA
AGTGGATTCACCTGTTTTTAAAAATAAACATTGATAAAAGGTG
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_012428.4
Summary:	This gene encodes a type I transmembrane protein belonging to the Ig superfamily. The protein is believed to be involved in cell-cell interactions or cell-substrate interactions. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2009]
Locus ID:	27020
MW:	41.8