

## Product datasheet for **SC212166**

### USP42 (NM\_032172) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** USP42 (NM\_032172) Human 3' UTR Clone  
**Vector:** pMirTarget (PS100062)  
**Symbol:** USP42  
**ACCN:** NM\_032172  
**Insert Size:** 1076 bp  
**Insert Sequence:** >SC212166 3'UTR clone of NM\_032172

The sequence shown below is from the reference sequence of NM\_032172. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CGTCTCTTTGAGTATGGCCAGGGTGATGAAAACCTCAGCCTCAAAACAAAAAATTCAGTAGTTATGATT
CAACGCGTTCAACAGAAGCCATCCCAGCCAGCTTAAATTATAAAGATAGACAATAACTCTGTTCCAA
TCTGCGTGGTGCTTCTTTAGTAAATACTGTACAGATTTTACCATGGAGAACTTTTTTTTTAGTTTTTAC
CTTTTCTTAATTACCCTTATTCCGAATGGACGAACACTTTCTACCACTGCTGACCATTGTAAAATACCG
TGATATAAAATCCCATTGAAATAATGCCCTGGAATAGAACATCTCAAATGCTGCTTAATTACAGACTCA
GGTCGATTACTTGTATTTTATGTAATGTTCCCAAGTTAGACATCTGGTGCAAGACCAACCGGGAGAC
CATGGAATTGTCAAAAGTACAACTGACAGTGTGTATTTAATTTAAAGACTTATTTAAAACTCACA
AGCTCTCACCTAGACTTTGGAGAGCAGTCTGTTTTCTGTAATGTCTGATACTAGAACTAATTTGCTTA
TTTTAGTTGTATTCAAGATTTGAAGATGATTTTTATAGACAAGTTCTGTTTTTGAACCTTTGTGGAAGT
TTCCAATCAATCAATTTCCAGTTATGATGAGTATTTACATTATGAATGTATAACCCAGACATGATTTG
TAAAGCCGACAGTATGTTTCTATTACACAACACTTTTTGATACAGCGTCTCTTGTCTTCACTGATACTG
GAGTCTCCGTTGTCTGCTTGGTCCCTTCGAGTTTCTAGTTACAGACACAATCATACTGTGATTTTATTT
TTAATATGGATATGCTATCAAAGTGTGATACACTTATAATTCAGTGGTCTGCATCAGGAGATGGAGTG
GGGAAAAGTGTATTTAATACAGTTTGTATCTGAATAATCTGTATGGTTTATACAGTTTGTGTTGTCAG
AGATGTTTAAAGTTTGTCTTTGTTTTTCTAAAGATTAAGAAAGCACTTGCCCCACTGTAATATACAG
CATGTAATTTCTATAGTATATAAATGGCAGCAAATCACA
AGCGGACCGACTTACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCC
CAACCTGCCATCACGAGATTTTCGATTCCACCGCCG
```

**Restriction Sites:** Sgfl-RsrII



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<u><a href="#">NM_032172.3</a></u>
<b>Summary:</b>	Deubiquitinating enzyme which may play an important role during spermatogenesis. [UniProtKB/Swiss-Prot Function]
<b>Locus ID:</b>	84132
<b>MW:</b>	41.5