

## Product datasheet for **SC212168**

### CaV1.3 (CACNA1D) (NM\_001128840) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CaV1.3 (CACNA1D) (NM_001128840) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CACNA1D
Synonyms:	CACH3; CACN4; CACNL1A2; Cav1.3; CCHL1A2; PASNA; SANDD
ACCN:	NM_001128840
Insert Size:	2000 bp



[View online »](#)

**Insert Sequence:** >SC212168 3'UTR clone of NM\_001128840  
 The sequence shown below is from the reference sequence of NM\_001128840. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GATGAAATGATATGCATCACCACCTTGTAGCCCCAGCGAGGGGCGAGACTGGCTCTGGCCTCAGGTGGG
GCGCAGGAGAGCCAGGGGAAAAGTGCCTCATAGTTAGGAAAAGTTTAGGCACTAGTTGGGAGTAATATTC
AATTAATTAGACTTTTGTATAAGAGATGTCATGCCTCAAGAAAGCCATAAACCTGGTAGGAACAGGTCC
CAAGCGGTTGAGCCTGGCAGAGTACCATGCGCTCGGCCAGCTGCAGGAAACAGCAGGCCCCGCCCTC
TCACAGAGGATGGGTGAGGAGGCCAGACCTGCCCTGCCCATTTGCCAGATGGGCACTGCTGTGGAGTC
TGCTTCTCCCATGTACCAGGGCACCAGGCCACCAACTGAAGGCATGGCGGGGGTGCAGGGGAAAG
TTAAAGGTGATGACGATCATCACCTGTGTCGTTACCTCAGCCATCGGTCTAGCATATCAGTCACTGG
GCCCAACATATCCATTTTAAACCTTTCCCCAAATACACTGCGTCCTGGTTCTGTTTAGCTGTTCT
GAAATACGGTGTGTAAGTAAAGTCAAGAACCCAGCTACCAGTGATTATTGCGAGGGCAATGGACCTCATA
AATAAGGTTTTCTGTGATGTGACGCCAGTTACATAAGAGAATATCACTCCGATGGTCGGTTTTCTGACT
GTCACGCTAAGGGCAACTGTAACTGGAATAATAATGCACTCGCAACCAGGTAAACTTAGATACACTAG
TTTGTTTAAAATTATAGATTTACTGTACATGACTTGTAATATACTATAATTTGTATTTGTAAGAGATG
GTCTATATTTTGAATTAAGTGTATTTGAACTGCAGCAATATCCATGGGTCCTAATAATTGTAGT
TCCCCACTAAAATCTAGAAAATTATAGTATTTTACTCGGGCTATCCAGAAGTAGAAGAAATAGAGCCA
ATTCTCATTTATTCAGCGAAAATCCTCTGGGGTAAAAATTTAAGTTTGAAGAAGTGCACACTACAGA
AATTTTTCTAAAATATTTTGTAGTCACTATAAACCTATCATCTTTCCACAAGATATACCAGATGACTATT
TGCACTCTTTTCTTTGGCAAGAGTCCATGATTTTGATCTGTACCTTTGGATCCACCATGGGTTGCA
ACTGTCTTTGGTTTTGTTTGTGTTGACTTGAACCACCCTCTGGTAAAGTAAGTGAATTACAGAGCAG
GTCCAGCTGGCTGCTCTGCCCTTGGGTATCCATAGTTACGGTTTTCTGTGGCCACCAGGGTGT
TTTTGCATCGCTGGTGCAGAAATGCATAGGTGGATGAGATATAGCTGCTTGTCTCTGGGGACTGGT
GGTGTCTGCTTAAGAAATAAGGGGTGCTGGGGACAGAGGAGCAACGTGGTGTCTATAGGATTGGAGTGT
CGGGGTCTGTACAAATCGTATTGTTGCCTTTACAAAACCTGCTGTACTGTATGTTCTCTTTGAGGGCTT
TTATATGCAATTGAATGAGGGCTGAAGTTTTTCATTAGAATGCACTCACACTCTGACTGTACGTCCTGAT
GAAAACCCACTTTTGATAATTAGAACCCTCAAGGCTTCATTTTCTGTCAACAGAATTAGGCCGACTGT
CAGGTTACCTTGGCAGGATTCCCTGCAATCAAAAAGATAGATGATAGGTAGCAATTTTGGTCCAAAAT
TTTTAATAGTATACAGACAACCTGTTAATTTTTTTTTTTTTTTTTTTTTTTTGTAAATAACAAACACCAC
TTTGTTATGAAGACCTTACAAACCTCTTCTAAGACATTTACTCTGATCCAGGCAAAAACACTTCAA
GGTTTGTAAATGACTCTTTCCTGACATAAATCCTTTTTTATAAAATGCAAAATGTTCTTTCAGAATAAA
ACTGTGTAATAATTTTATACTTGGGAGTGCTCCTTGCACAGAGCTGTCATTTGCCAGTGAGAGCCCTC
ACGCGT AAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_001128840.3](#)

**Summary:** Voltage-dependent calcium channels mediate the entry of calcium ions into excitable cells, and are also involved in a variety of calcium-dependent processes, including muscle contraction, hormone or neurotransmitter release, and gene expression. Calcium channels are multisubunit complexes composed of alpha-1, beta, alpha-2/delta, and gamma subunits. The channel activity is directed by the pore-forming alpha-1 subunit, whereas the others act as auxiliary subunits regulating this activity. The distinctive properties of the calcium channel types are related primarily to the expression of a variety of alpha-1 isoforms, namely alpha-1A, B, C, D, E, and S. This gene encodes the alpha-1D subunit. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012]

**Locus ID:** 776

**MW:** 75.4