

Product datasheet for **SC212170**

CCDC102B (NM_001093729) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC102B (NM_001093729) Human 3' UTR Clone
Symbol:	CCDC102B
Synonyms:	ACY1L; AN; C18orf14; HsT1731
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001093729
Insert Size:	1082 bp
Insert Sequence:	>SC212170 3'UTR clone of NM_001093729 The sequence shown below is from the reference sequence of NM_001093729. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACTGAACCTCAGGCACTTGCAAACTGGTATTTTTTCACAAAATATGCTGAATTAAGATTAGGGCCTT
AAAGACATTTCCATATCCTTTTCTTAAATATCAGTAAAATTGTTTTTAACTAGAAATATTAATGAA
AAAAACGTAGACAATACACAAATTAATGGGCTTCTTCACTTCTTCTAATTTTGCCTAACAGATACTGC
ATATTCTCAAAAAGACAATTTAAATGTCATTTAAAAACAACCTTTAATTCTAAGATGTGTAATATTTTG
AAAGTCAAAAAGGCTTTTCAGAATACTTTTACATAAAATCTGAAAGAGTTATAATATCGGTAAGAAAA
AGTAAGTTGAAAACCATACAAGACGCTGGGTCTTAATAAGAAAACCATTGACTTTAAGTATAAAGTAC
TGGTTTGTAAATAATTGGTAACTTTTATGTACGTGTTGTCTATGTGGTGGGGATGGCAGGTTGTAT
TAACAAAAATGAATCATTCTAGAGGTGTAACAATACATTTCTTATATAATTTTATAAGTCATTTCTAAT
CTTTGTATAAACAGAAGTGAGCAGATGAATCAGAAAAAGTGTTTTGTATTTTAAAGTAACAGATAAC
CAGTGATTGAATCTAAGACAGGCTGTAAGCATCGCTGAGAACTAAAAGGACTTTTGACTTTTATCTGG
ATAGACATTTCTACAGTAAAATCATGGAAGGCATCAGCATTGCAAAGTAGCATCTAGGTAGAAATCAG
GCCAAATTAAGCTGTGGTTTCCCTCTGAGTAGTGGAATAGAGAAAATTAGGAAATTGTGGTTATGTG
AATATTTCTTTAAACTTTTATGTACATTATAGTTTATTGCTTCATATTTAAGTTTAAAGTTTAAAGGTA
AAATGTTATTTTGAACAAAAAGACACTTATAATTTTCCATACCTATTTTCAACTGAAGGCAACTGTAA
GATTTAACTCAGTCAATAACATACTGGTTTTACTCATCTCCCCTCCATTGATTAGCCAAAAAAAATG
AAATCTTACTAATTCATTATTGAATAAAGACCACTTTTATCAGAGTC
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001093729.2
Locus ID:	79839
MW:	42.1