

Product datasheet for SC212182

WTAP (NM_152857) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	WTAP (NM_152857) Human 3' UTR Clone
Symbol:	WTAP
Synonyms:	Mum2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152857
Insert Size:	1076 bp
Insert Sequence:	>SC212182 3'UTR clone of NM_152857 The sequence shown below is from the reference sequence of NM_152857. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AGTGCCTGGAAGTTTACGCCTGATAGGTAACAAATCATACTCCCCAGTCAAGACTTCCCTGACAGTCC
CACTACGAGAAAGCTGTGGTGGGACAGCAAGTACTCGTTTCCACACCAAGACTCAGACTTTTTGAGCC
AAAAAAAAAGCCACATTCTTACACTGTCCAGCTTGAATGGTTAATGTAAACTTACCAGATGAACCTTG
TGTTTCAGCTTTTTCTTTCCCTTCCCTTGCTTCAGAGGCCTGATGGCGTCGGACTATTCGAAGA
AGTGGCCACCTCCGAAAAATCCCTTCTAGAACATGTAGACACTTGAGAAATGTTTCTGTTGAAGAA
AATAGAGGGAGAAACAGAAGTCTTAAGTCTGTGGCACACTGTGTCTTCAGACAGTTTGAAGGAATGAAA
ACCTAGAGATTTTAAATCATGAATTGAACATGTAAATTCAGTAAATGTAAAAACGGAATATGCATC
GCTCTTAACCTTGAGCATAGTGACTTAGAGACACTGTGTATCAGTTTTGCCAATAAGACTGTGGACTTC
ATGATTGTTGTTGAACCTCTGGGTCAAACTCAAATGAGGTGAATTTGCCTTTAAAGGGTTTATTTGC
TGAGAACCAACTTTCAATAGTCATGAGAGAATCAAATAATAGATGTCCGTACAAGTAGCGCATATATTT
AACCATTTAGTTTGGGGCTCTATATTACTTGCTTGAGCCTTAATCAATGTGGTTTTATTCAATGGTTTG
TTCTTTGAATGGTTGCAAAACTGTAGATAATCTTACTGAGGACTGTACAAACATGAAGGTGTGGTATC
AAACTTCAGGTTGAAACTGTTTGAAGCATTATAAACATTCAATTCACAAGTATGTTATGAAGGATATT
AGCTGTGATGAGACTCACTGCATTATTTTTTTAGTGAATTTTATGAAATCCCCGTTCATTCAACAGG
CACATGTTTAAAGAGCTTTGTCGTTGGTGTAAATGGGGGAATGTGTTCTTCATTGATTTGGGCCTT
TTGTATTGCACTCTTGATATTAAATTAATGTGCCTTGAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152857.3</u>
Summary:	The Wilms tumor suppressor gene WT1 appears to play a role in both transcriptional and posttranscriptional regulation of certain cellular genes. This gene encodes a WT1-associating protein, which is a ubiquitously expressed nuclear protein. Like WT1 protein, this protein is localized throughout the nucleoplasm as well as in speckles and partially colocalizes with splicing factors. Alternative splicing of this gene results in several transcript variants encoding three different isoforms. [provided by RefSeq, Jul 2012]
Locus ID:	9589
MW:	41.4