

Product datasheet for SC212198

C7orf57 (NM_001100159) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C7orf57 (NM_001100159) Human 3' UTR Clone
Symbol:	C7orf57
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001100159
Insert Size:	1055 bp
Insert Sequence:	>SC212198 3'UTR clone of NM_001100159 The sequence shown below is from the reference sequence of NM_001100159. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TCTGCATCAACACCAGCAGAGCTCAAATAAATCCTGATGCAATATGTATTTAGGATAATTTTTAAATGG
CTAAATATGACATGACTGTATTATAGCTATATTTCTGAGGCTTTTTTGTATTTTATTAATATAGACTC
TCATTGAATAATTTACCTTGCAGCAGATTTGACATTGCTGCATAGAAGGGGCAGGTGTTGGTTTTGT
TCTTGCATTTCTCTGGGATGTGAAACACTTGGCTAACAAACACAATAAGGCCTCTGGTACCCTCTGCC
CTGCTGGCCTGGTGATGTCCCAGTGTCTCAGAACTGATCACTTGCTCCTCTATTGTGTGAACCACTCT
TAGCAGGCATTGCTACTGGATTCTTAAGTAGAGAAAACCTTACAGAACAGGTGTTAAGCCTTTTCCAATA
TTAGTATTGTAGTCTTCATAGAAATCATTTCTAGTTACAGAAAAATTTAACTCCAGAATATGGTGATA
CAGTTTATGTTAAGAAACATATTTTAAATTTATTATGCCATGAATTTATAAACATTTAATATAAACATT
ATAACATTTGCTGTTTTAAATTTTGTCTATTTTTTACAAACTGATATTCAAAATATTAACATTTGT
TATTTTGTATAAACCTTTAAGTAAAAATATGCTGTTTTAATTTAATTTAGTAGAATTAACAAAACT
TTTTAATAGTCAAACTGAATTACCTTGTTTTATAAATGTTAAATGTTTATTTAAGTAATTAGCTTA
AAAATATAGTTTTATGTTTTGTACATTAATTGATTAATAAAATAAATTTTATGTGAGAATGTTTTCA
AAAAGCCAAATAAAGTCACAGATATTTTAGCTTCTTTCCAGAATAAACTTATCATGCATTGCTTACC
TAAATTTAAGTTTCAAATTAGTTCTCTATATGAGATACTTAATCGAGCATATATATTTATATTCATTG
ATCAAACTGATATGCTACAGGAGAATGTATGCAATGTGTAATAAATGCAATGGCTGTATGCATTGAAT
AAAAATAATATACAATTAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001100159.3</u>
Locus ID:	136288
MW:	41.6