

Product datasheet for SC212202

PAXBP1 (NM_016631) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PAXBP1 (NM_016631) Human 3' UTR Clone
Symbol:	PAXBP1
Synonyms:	BM020; C21orf66; FSAP105; GCFC; GCFC1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_016631
Insert Size:	1081 bp
Insert Sequence:	>SC212202 3'UTR clone of NM_016631 The sequence shown below is from the reference sequence of NM_016631. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAATTTAAGTCTTTGATCGAAGGAAAAATAGTCTATGAGACTGAAAAACCTATTTGCTTGATTACCATT
CCCAATGTGTAATTTTGTATATATGTAAAGTTTCTTAACTGTAAAATAGCAATTCTTGTTTAAATAC
AAAAAACATTATATTCAATACTGTGCTTAAATGATATTTCTTCAGAATGAAAGGAATTCAGATTATT
TCTAGGAGTCTATGAATGTTTTCTTACTGCCTGTCATACTTTGTTTACTGAAGTATTTGTATGGTCA
TTTAAATTAACCCCTCTCAGTTAATTGTCCCTGTAAACGATGTGTGCAGTGTAAATTGTGTAATTATGC
ATATATATATTTATACCACCATGGACTTACTTTACACTGGAAATTTGTGGTGTTTTCCACAGAGGAA
GCCTTTTAAACAGCAGAAAATTATCCGGTAGCAATTTGTCTTAGGGAAATTAATACTGTTGTAACAAA
AGGGCTGGCAAGTGACTAAATGTTACCTTTTGGTTAGCGTTGCTATCACTGTAGCTGTGCTGCCT
TAACATTGAGAGTTGAGTTGTGGACTTAGTCATTGAGTGAATGTTGAACCTTTCTTCTTGAATCAGCGT
TAGAGTCAAAATGTGGGTTTAATGTGAATGCTAAATGAATATGCAACATAGTTTCAATGGAACATAGAA
ATTTCACTTTTGAATAATAAACCACAACTAATTGCATGGTTTTGGGCTGACCGTCTGTCTCATT
CATGATTTTACTCCCTGCCAGCTTGCAGATACATGCAATTTGAATTTAGTAGTCTGTCATATTTAAAA
TCCACATATGTAGTTAAATATGCAATAAAGAATACGCCTTGCTGATTAGCCACACCCAGGCAATTTA
ATATCCATGGTCTTCTATAATGTTTATACCACATTGTAAATATTAGAGAACTTGTTTTAAAGTTTATG
CTTTCTGATTCTGTGATTTTCAAAGAGGAATTTTTTGGAAATTAAGGACTATCGTTTGATAATTATA
TATGTAATTAAGATTTAGCCAAGCATTAAATTTTAAATATA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_016631.4
Summary:	This gene encodes a protein that may bind to GC-rich DNA sequences, which suggests its involvement in the regulation of transcription. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2009]
Locus ID:	94104
MW:	42.2