

Product datasheet for SC212204

FOXD2 (NM_004474) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FOXD2 (NM_004474) Human 3' UTR Clone
Symbol:	FOXD2
Synonyms:	FKHL17; FREAC-9; FREAC9
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004474
Insert Size:	1098 bp
Insert Sequence:	>SC212204 3'UTR clone of NM_004474 The sequence shown below is from the reference sequence of NM_004474. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GTCGCCGGCCTTAGTGGCTGCCACTTCTACCGCGAGCAGGCCAGGGCCGGTTAGGTCCGCACTCCTCA
GCCTCTCCCGGGAGTTCCTGCGGTCCCAGCGAACTCAGGGAGTCTATTTATGAAGTCTCCAGACCTTG
GGCCGGCACGCGTGACACGGCACTTCAGGCTCCACGCACAGAATCTCGCAGATAGTTGGGACTAAGCGG
GCTCTATCGCTCAGGGCGACAGGCCCGGGGCTACGCGAAGAAGTCGAGGCCAAGATTCTTTACAGTTT
GAGAAATAAAAGCAGGGGGTGGGGGCTTCGTTTTTTTCCCTGCCCTCTGCGCCTCTCGGGGAACACATT
CCGGGAGAGATGCCTGGCCAGGCTCCACGGATCCCGCCAGAAACACCAACAGAGGGTCTCCCTTTCTGC
CTTTCCCTCTCACTTCTTCCCCAACGTTGAGACCCTGCTTGTCATATTATAATTTAAAGACATCTA
TTATCTGCTTTGTGCTTAAAGAAAAATTCAACCTTTTTTTTTTTTTTTTTTTTGTGTTCTCCAA
GGAAGTTCGTTTCTCTGAAGCCTAAACCAGTGTCTACGCAGGCGGAGCTGAACGGAGAGGTGAAGCAG
GGGTCTTTATATTCCCTGCAGAAACCCTGGATCCCACTCCACTGAGGAAATCCTTTCTTTTGTGT
TTGGGTGCTTTTAAAGGAATTTCTTTCTTTCCCTTGGCTGGCTGGGTCCAGCCATCTCCGGCCAGCAA
CCCTGTCAAGTGTCCACCCTGTCCAGGCCCTGTGATATTTATTGTGGCTAAAGCAATGAGATTGGGT
GGAAGAGAAGGGTCTGGGCCAGGGACTGGAGAGGTGGGGGTGGAGTGGGCAAGGGCACACCCTGTGGCA
GGGGGTCTATAATTTTGTGTAATATCTGTACAGCGCGCTGCTCTCGGCGTTCTGTGCACTGGGTG
TGCATTGATTTAACTTCTCCACTGTGGGAAGCTCTGTGAAATCGTCTGATGGTCTGTGTGAAGAAA
AGAAAAATCTGTCTGCCGCTATTTCCATTGGTGTGATTGAATAAATGTATGTGAACCTCA
AGCGGACCGACTTACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAATGACCGACCAAGCGACGCC
CAACCTGCCATCACGAGATTTGATTCCACGCCGC

```

Restriction Sites: SgfI-RsrII



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_004474.4</u>
Summary:	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined. [provided by RefSeq, Jul 2008]
Locus ID:	2306
MW:	39.3