

Product datasheet for SC212212

Prokineticin 2 (PROK2) (NM_001126128) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Prokineticin 2 (PROK2) (NM_001126128) Human 3' UTR Clone
Symbol: Prokineticin 2
Synonyms: BV8; HH4; KAL4; MIT1; PK2
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_001126128
Insert Size: 1099 bp
Insert Sequence: >SC212212 3'UTR clone of NM_001126128
 The sequence shown below is from the reference sequence of NM_001126128. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AACCGATTATTTGTTTAGCCAAAAGTATCGCTCTGGAGTAGAAACCAATGTGAATAGCCACATCT
TACCTGTAAGTCTTACTTGTGATTGTGCCAAACAAAAATGTGCCAGAAAGAAATGCTCTTGCTTCTCT
CAACTTTCCAAGTAACATTTTTATCTTTGATTTGTAATGATTTTTTTTTTTTTTTTATCGAAAGAGAAT
TTTACTTTTGGATAGAAATATGAAGTGAAGGCATTATGGAAGTGTCTTATTTCCCTGTTTGTGTTT
TGGTTTGATTTGGCTTTTTCTTAAATGTCAAAACGTACCCATTTTCACAAAAATGAGGAAAAATAAGA
ATTTGATATTTTGTAGAAAACTTTTTTTTTTTTTTCTCACCACCCCAAGCCCCATTTGTGCCCTGCC
ACACAAATACACCTACAGCTTTTGGTCCCTTGCCCTCTCCACCTCAAAGAATTTCAAGGCTCTTACCTT
ACTTTATTTTGTCCATTTCTCTCCCTCCTCTTGCAATTTAAAGTGGAGGGTTTGTCTCTTTGAGTTT
GATGGCAGAATCACTGATGGGAATCCAGCTTTTGTGTCATTTAAATAGTGAAGAGTGTATATGTG
AACTTGACACTCCAACTCCTGTCATGGCACGGAAGCTAGGAGTGTCTGTCGACCTTCCAAACCTGT
CACTCAAGAGGACTTCAGCTCTGCTGTTGGGCTGGTGTGTGGACAGAAGGAATGGAAAGCCAAATTAAT
TTAGCCAGATTTCTAGGTTTGGGTTTTCTAAAAATAAAGATTACATTTACTTTTACTTTTTAT
AAAGTTTTTTTCTTAGTCTCCTACTTAGAGATATTCTAGAAAATGTCACTTGAAGAGGAAGTATTTA
TTTTAATCTGGCACAACACTAATTACCATTTTTAAAGCAGTATTAAGTTGAATTTAAACCTTGTGTTG
AACTGAAAGGTCGATTGTAATGGATTGCCGTTTGTACCTGTATCAGTATTGCTGTGTAATAAATCTGTA
TCAGAATAATAACAGTACTGTATATCATTTGATTTATTTAATATTATATCCTTATTTTTGTCA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001126128.2
Summary:	This gene encodes a protein expressed in the suprachiasmatic nucleus (SCN) circadian clock that may function as the output component of the circadian clock. The secreted form of the encoded protein may also serve as a chemoattractant for neuronal precursor cells in the olfactory bulb. Proteins from other vertebrates which are similar to this gene product were isolated based on homology to snake venom and secretions from frog skin, and have been shown to have diverse functions. Mutations in this gene are associated with Kallmann syndrome 4. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	60675
MW:	42.8