

Product datasheet for SC212214

C15orf39 (NM_015492) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C15orf39 (NM_015492) Human 3' UTR Clone
Symbol:	C15orf39
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015492
Insert Size:	1083 bp
Insert Sequence:	>SC212214 3'UTR clone of NM_015492 The sequence shown below is from the reference sequence of NM_015492. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CTGCTGGACAGCCAGAGCCATCACCTGTAGCACTGGTTGCCAGTGCTGTGTATAGCAGTCACTCTCC
ACCTTCCCTTCTGCCTGCCAGCTGCCCCGGGGCCACGAGTGGATGCTGGGGCTGTGGCTGCTCCCT
GGAGGGTTCCATCTCTGACCTGTGGCCATTGAGGTGGGCTGAAGAGCCCTGAGCTTTTAAACGTG
AGGGTCTTTATTGGATAGGACTACTCCCTATTTCTTGCCTAGAGAACACACATGGGCTTTGGAGCCCGA
CAGACCTGGGCTTGAATCCCGGCTCGTGTCTTCTGCTGCAGGACCTGGGCAAGAACTTCACCTCTGCTG
AGCCCTCATTCCCATGTGTAAATGGGACAACGCAACCTACCTCACAGGTTGTTGTGGGGATGCTGC
CTGATACATACCCTGTCACCATTTGGTCTCTGCTTCTCTCTGGGACAGGGCCTAGAATTGGAGGCAGA
GAACCTTCCTATAGAAAGTCTTCGTGTGTCTAGGACTTGGCTATCGTAGAGTGGTACCTTAGGCAGTG
GATGTGACTCACACTTTCAGGAGTCACCCCCAGCATTGGGGTTGGGTGGCCCTACTCCAGCCTGGA
GCTCCCTGAGGGAGCCTGCACTCCCTGCTCCCAATCCCGCTACTGGTGCAGGGATGCAGCCTGGAGCT
GGCGTCTTGTCTGGGCTGCTGCTGCCGCCACCCAGGAGGCCAGGCCTGCTCTGAATTGACATC
AGTGCTTCCCTGAACTGCCTCCCCACCCCTGGCATTATCCAGGAACTTATGTTTTCTAGAAGCTAA
GCAGCTGCTGGGACTCAGGGACTGGTGCAGGTAGGCTGAGTGGCAGCTCAGTCCTAGAAGGTCTCTGAA
GATCTGGACTGAGGACCCTGCTGCTCCCAAGCCAGAGCCCATCAGCCAGGCCTGCTGTGAGCCACCTG
CCTGTGGAGTGCTGAGCTCAACCAAAGGCTGGCAAGCTCTGGGCCTCATTTAAGGGATTCTGATGAGCC
GATGGGCCTGGAGGCAGCCATTAAAGCATCTGGCTCGTTTTTGGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015492.5</u>
Locus ID:	56905
MW:	38