

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001130864.2
Summary:	H2A.Z-specific chromatin binding protein which may play an important role in the neural crest stem cell migration and differentiation during early development. Also required for proper mitosis progression.[UniProtKB/Swiss-Prot Function]
Locus ID:	114825
MW:	40.6