

Product datasheet for SC212222

SLC66A3 (NM_152391) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC66A3 (NM_152391) Human 3' UTR Clone
Symbol:	SLC66A3
Synonyms:	C2orf22; PQLC3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152391
Insert Size:	1058 bp
Insert Sequence:	>SC212222 3'UTR clone of NM_152391 The sequence shown below is from the reference sequence of NM_152391. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TACCGGAAGACCGCTATAAAGGCTGAATGATGGATACATTATTCCTTCACACAGTGGATTTTGAGTAAC
TGAACCAAAGGAAAAAGAAGCTCTTTGCTAAATTAAGGTCTTTTATAAATTTAGTAAATCAGTTTATAA
TCTTTAAAGCCAAAGGTTTTTTAGACTTGAAAGAAAGAGCCACTTAAATTCCTGTTTAAAAATACCAA
TTTGCCTCCTCCTCCTCACTTCGTTAGGTTATGGTAGTGCTCAGACATCTGCAGTGTTGAGGCCAGTC
ACTGTTGGAAGTCATCCAAGAAGCCATTTTGAGGCCATTTTGAGCCTTACTCTTAAGTTCTCTATGAA
GAACTACATTGATTTGTTGGCTTTCAGAATCTTTTAGGAAATAAATCCTCTCCAGGACAAAAATGAACA
TGAATGGAGTGGCATTGTTTCCAAGTCAGAGGTGGGCACCTATAATAATGACTAGGGTTCACCTTCT
GGGACTGATGTTTAATTGTAACACAGATACAACAGGGTGGCCTTGTTGTGTATAATACGGTATTATACC
CGCATGTGCTCTAGCAAGGATACCAAGGCAAGCATACATGTAGCTGGCTTGAGTTTGTACCAAAACAGT
CCTTCACTTTGCACTGTGCCTTAAGTAATTACTAACAAAAGTACTAGGATTAGCTGCAATCTCTACT
TTCGATGAGGAAATCCCAGTAAGCTTTCTGATTCAAGTACAATGCTGCCATTTTTTAAAGGGCCACAAC
TATAGAATTACCACTGTTGGAATTTGGTACAAAATATGTTTTGTCTATTGAAAACATACACGGTAAATG
GTGTTGTTAGGTAGGTTCTGTCCAGTTCTTAGGGACTTTTTTTCACATTATAGCATTTTTACCCTAAACA
TGATGTTGAGATTATTATATACTGTATTTTCTTCTAAATTAACCTAATGTTTAAAACTCACTTTCCC
CCTTTAATTGAAGGCATTGTTTTGTAGATGCAGTAATGATGTTTACCAGAGATTATTGTTTCCTATGC
AAAATAAATTTTCATATTTTGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_152391.5</u>
Locus ID:	130814
MW:	40.5