

Product datasheet for **SC212258**

PTGES3L (NM_001142654) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PTGES3L (NM_001142654) Human 3' UTR Clone
Symbol:	PTGES3L
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142654
Insert Size:	1091 bp
Insert Sequence:	>SC212258 3' UTR clone of NM_001142654 The sequence shown below is from the reference sequence of NM_001142654. The complete sequence of this clone may contain minor differences, such as SNPs. Red =Cloning site Blue =Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**CGCATCGC**

TCTGACAGTGCTGATGATGCAACAAGTAAT**TA**CTTTCTGTGACGCAAAGCTGGGAAGGCAGCTGTGGCT
 ATTTCCAGTTGTTCTAGAAAGCTAGCGCCTTAGGCCTTTGTGACGCTTTGGCTGTGCACCAAGTTCTTC
 TGAGATCTCTGAACCTTCTAGAAAGCTCTTTATTAAGGGTCGTCCTTCAAGCTTGTCTTCCTCCTAGGCA
 CTTGACTATGGCTCCACATCAGATGCAGGGTCGGGGAGGGGTAGACTGGTAGGGCTTCTGATTTCTTACT
 GATTGGTTTCTAATCTCTGTGTTTAGGCCACTTGATTATTTGAGTATGTGTGTATGTGGCAAAACAAA
 AGCCAAAGAAGATGCAGATTATAGAAATCTCTACTATGTGTGCATATGGGGCAAAGCAAAAGCCAAAG
 AAGAAACAGATTATAATAAATCTCTAGGAAGTGTCTGTCTGTTTCTTCTTGTACTTATTCCTCCTG
 GTGGCAAGCCTCTTCTCCCTCTAGTAACCAAGCTGGGGTGCTTCTCTTGCCAGATTTGTGCTCTGATT
 TCTTACAGCTAATTCTGTCATCTCTAGAACCAGCCACTTCAGCTTCTGCACACTCTAGAAATCTCTTTA
 TACTTTTTTTTTTTTTTCCGAGATGGAGTCTTGCTCTGTGCGCCAGGCTGGAGTGCAGTGGTGCAGTCT
 CTGTTCACTGCAACCTCCACCTTCTGGGTTCAAGCAATTCTCCTGCCTCAGCCTCCCAAGTAGATGGGAC
 TACAGATGTGCACCACCAGGCCAGCTAATTTTGTATTTTAGCAGAGATGGGGTTTACCATGTTGGC
 CAGGCTGGTCTCGAACTCCTGACCTCAAGTAATCCACCACCTCAGCCTCCCAAAGTGCTAGGATTACAT
 GTATGAGCCACCACCCGGCCCTCTTTATACTTTTAAACTGATTGTTACTCTAGTTAATGTCTATGTT
 TATTTGTCCAACAGCTATTCTATTTGTCCATAAATTTTTGAGTGCCTAATATGTACAGTTGACTGTAA
 CTTAAACAAAATAATGCATTGCCTATAGACTTTGCTCTTG

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001142654.1</u>
Locus ID:	100885848