

Product datasheet for SC212281

GLB1L2 (NM_138342) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GLB1L2 (NM_138342) Human 3' UTR Clone
Symbol:	GLB1L2
Synonyms:	MST114; MSTP114
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_138342
Insert Size:	1086 bp
Insert Sequence:	>SC212281 3'UTR clone of NM_138342 The sequence shown below is from the reference sequence of NM_138342. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CACCTGGGCAGGAACAGTACATTAAGTACGGTGGCACCCTCCTGCTGGTGCCAGTGGGAGACTG
CCGCTCTCTTGACCTGAAGCCTGGTGGCTGCTGCCCCACCCCTCACTGCAAAAGCATCTCCTTAAGT
AGCAACCTCAGGACTGGGGCTACAGTCTGCCCTGTCTCAGCTCAAAACCTAAGCCTGCAGGGAAA
GGTGGGATGGCTCTGGGCTGGCTTTGTTGATGATGGCTTTCCTACAGCCCTGCTCTTGCCGAGGCT
GTCGGGCTGTCTTAGGGTGGGAGCAGCTAATCAGATCGCCAGCCTTTGGCCCTCAGAAAAAGTGCTG
AAACGTGCCCTTGCACTGGACGTACAGCCCTGCGAGCATCTGCTGGACTCAGGCGTGCTTTTGCTGG
TTCCTGGGAGGCTTGGCCACATCCCTCATGGCCCCATTTATCCCCGAAATCCTGGGTGTGTACACAGT
GTAGAGGGTGGGAAGGGGTGTCTACCTGAGCTGACTTTGTTCTTCTTCAACCTTCTGAGCCTTC
TTTGGGATTCTGGAAGGAACCGCGTGAGAAACATGTGACTTCCCTTTCCCTTCCCACTCGCTGCTT
CCCACAGGTGACAGGCTGGGCTGGAGAAACAGAAATCCTCACCTGCGTCTTCCCAAGTTAGCAGGTG
TCTCTGGTGTTCAGTGAGGAGGACATGTGAGTCTGGCAGAAGCCATGGCCCATGTCTGCACATCCAGG
GAGGAGGACAGAAGGCCAGCTCACATGTGAGTCTGGCAGAAGCCATGGCCCATGTCTGCACATCCAG
GGAGGAGGACAGAAGGCCAGCTCACATGTGAGTCTGGCAGAAGCCATGGCCCATGTCTGCACATCCA
GGGAGGAGGACAGAAGGCCAGCTCAGTGGCCCCCGCCCCACCCCCACGCCGAACAGCAGGGGCA
GAGCAGCCCTCCTCGAAGTGTGTCCAAGTCCGATTTGAGCCTTGTCTGGGGCCAGCCCAACACCT
GGCTTGGGCTCACTGTCTGAGTTGCAGTAAAGCTATAACCTTGAATCACA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_138342.4</u>
Locus ID:	89944
MW:	38.3