

Product datasheet for **SC212288**

P2Y2 (P2RY2) (NM_002564) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	P2Y2 (P2RY2) (NM_002564) Human 3' UTR Clone
Symbol:	P2Y2
Synonyms:	HP2U; P2RU1; P2U; P2U1; P2UR; P2Y2; P2Y2R
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_002564
Insert Size:	2000 bp



[View online »](#)

Insert Sequence: >SC212288 3'UTR clone of NM_002564
The sequence shown below is from the reference sequence of NM_002564. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AGCGAGAACACTAAGGACATTCGGCTGTAGGAGCAGAACACTTCAGCCTGTGCAGGTTTATATTGGGAA
GCTGTAGAGGACCAGGACTTGTGCAGACGCCACAGTCTCCCCAGATATGGACCATCAGTGACTCATGCT
GGATGACCCCATGCTCCGTCATTTGACAGGGGCTCAGGATATTCAGTCTGTGGTCCAGAGTCAACTGTT
CCCATAACCCCTAGTCATCGTTTGTGTATAAGTTGGGGGAATTAAGTTTCAAGAAAGGCAAGAGCTC
AAGGTCAATGACACCCCTGGCCTGACTCCCATGCAAGTAGCTGGTGTACTGCCAAGGTACCTAGGTTG
GAGTCCAGCCTAATCAAGTCAAATGGAGAAACAGGCCAGAGAGGAAGGTGGCTTACCAAGATCACATA
CCAGAGTCTGGAGCTGAGCTACCTGGGGTGGGGGCCAAGTCACAGGTTGGCCAGAAAACCTGGTAAGT
AATGAGGGCTGAGTTTGCACAGTGGTCTGGAATGACTGGGTGCCACGGTGGACTTAGCTCTGAGGAGT
ACCCCCAGCCCAAGAGATGAACATCTGGGGACTAATATCATAGACCCATCTGGAGGCTCCATGGGCTA
GGAGCCAGTGTGAGGCTGTAACCTATACTAAAGGTTGTGTTGCCTGCTGAGCTGTGCCCTATTGTGTGG
TCGGGGGATGAGGATATGGCAGGGAAGCTTTCACCAGCCACACAAGGTCCTTTCTCAATCCGTTCCC
TTCTGCCACCTGCCTTCTCACTAGCTGTCTCAGGAGTAGTCTCATATCAGGGATCCTCTCTCCAGGCCC
CAGGTCATCCCCACTGTAAAGCCAGTTGGCTTCTGTGCCTGACTCTGTGCTGAGCACAGAGAAAAGTC
AGGTGCAGTCCCAGTCTTGAGATTTCCAGTTTAGGTGATGGCCAGTCATGTGCTGGTAAATGTTAAG
CCATTTAACTGGAGCTCCGATTTAACTGGGAACCCCTTCTTTGTAGAAGTGGCCATTTCCATGGTGCAA
AACTCTTATGGCCAAATTCAAACATAAACATGATGTCAACTTCCTGTAAAATTCCTGAAAATTTAA
CAATTGGTTCTTGCAAGCAGGTACAAGTCACTCTAGCACACCACTGGATAATGCCGAGTGGCTGGGGCT
GTGAGCCAGGGGTCAGGGAAGGTTTCTGGAGTGGGAGACTCCTGAGCTGTAGCTTCTGAGAAAAGCA
GCTCACCGGAAGAACATCCACACACAGGATGCATCCAGGCCAAGACATGGGGCAAGAGGCGAGGTGT
GCATGAGCCATAAGCAGGCTGGACAGCTGACTCCAGGCTCAAGGCAGGGTTGGGGCTGGGTACAAAGGA
GGCCAAGTTAGGGCTCCCTCCTTGCCCTGACCCTGAGGCTTCTTTGCACTGGGGGTAATATGAAAGA
GATTCACCTCCTCCTCTGTCATCGTCTGCCTTCTGGGAAAATCCAGAATGGCAGGGTGGTGTCTCAGG
CTGGGTGAGGACTTAGTATGGGGGAGATGGGTCTCACCTATGATAGGTTCTGAGTCTAGCCAGGACCAC
TCTGGGCTGACGTGTGGCGCTCAGCTGGACAGGGCCCCGCCCTAGCCTTTGGAAAGGACAGGGCAAAG
CTGACAGGCTCACTCTTGATCTCAAGGACAGGACTCCCTCCTCTCCCTCTGGGAGAGCCCTCGCCCT
GTGGCCCACTCTGCCTGCCCCCTCTGACCTCTCCACCCTTCCCACTCTGCCTTGTGAAAGGCAGGACAT
CCCCAAAGCTTGCTGTATTTGGAGCCTGACTCCACAGCGCCCTCATGTGACAGAGACCCATCACAGACC
ATGACCCAAATGATTACTCTGTACTGTGCCAGGTGGGTGACCTGCCCAGAATAGTGTGAGCTATTACC
TCTCACCTTCTAGGTTCTATACTTCTGTTAATGTAGCCGATGCCTTCTGAGTTTTTTTTTCTTTTT
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_002564.4](#)

Summary:	The product of this gene belongs to the family of P2 receptors, which is activated by extracellular nucleotides and subdivided into P2X ligand-gated ion channels and P2Y G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor, found on many cell types, is activated by ATP and UTP and is reported to be overexpressed on some cancer cell types. It is involved in many cellular functions, such as proliferation, apoptosis and inflammation. Three transcript variants encoding the same protein have been identified for this gene. [provided by RefSeq, Mar 2013]
Locus ID:	5029
MW:	72.7