

Product datasheet for SC212295

GTP cyclohydrolase 1 (GCH1) (NM_001024024) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GTP cyclohydrolase 1 (GCH1) (NM_001024024) Human 3' UTR Clone
Symbol:	GTP cyclohydrolase 1
Synonyms:	DYT5; DYT5a; DYT14; GCH; GTP-CH-1; GTPCH1; HPABH4B
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001024024
Insert Size:	1084 bp
Insert Sequence:	>SC212295 3'UTR clone of NM_001024024 The sequence shown below is from the reference sequence of NM_001024024. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAAGAGTTCCTGACTCTCATTAGGAGCTAGCTTCATTCAAGTGTGTGTCTGCAGAACCATAGCTTCCAC
GCACCTGAACGAGCACAGAATGAACTGACGGTGGAAAGACATTATGAGCTGTGTCCAACGTTTTAACCAA
AGCGTATCGTACCAACGATCTGTGAAAAATGCACTGGAAGCTTCTGGTCCCGGTTTCTTTGTGGTCTAT
GTGGGTCTGTCTCATTGTAACCTCCGTATAGATGGTATAGGTATTTAATCCTGGAAGCTGTTGCCTT
ATTAATGATTATCTTAAATTTCTCCATTGGGGCAGCGTGGGCCAAATTAACAAACAAACAAACCGCAA
CTCCTCCACAGAAACACAAACACAGTTATTCCATGAAGTTTAGTATTTGGTTGACATAGTCTCTCAA
ATTCATCCCATACCTAAAAGTAATAACTTTGATGCTTGCTTTAACTTTAGTCCCATCTCTGCCACTT
TGATGCTATTTGGGTTATGATGGGGCAAGATGGCAGAGGTATTGGGTTTTTTGTTTTTTCCATTCCCT
CTCTACTTCTGTTTCTAGCTTTTTCTTCTGGAGTTAAGTACAGTGATGGTTGGCTTGAGTACTTTT
TTAAATCTAGCCAGTATAAACATTAGCCTGCTTAATATTTAGACATTTATAGGTAGAATTCTGAGCAC
TCAACTCATGTTTGGCATTTTAAAGTAAAAACAAGTGTGACTTCGAGGACCAAGAAATTTGTCAGCTAT
ACATTTATCTTTATGAACCTATTTATTTCTTTTTAATGACTCGTTGTTCTAACATTTCTAGAAGTG
TTCTTATAAAGGTCTAATGTATCCACAGGCTGTTGTCTTATTAGTAAATGCAAAGTAATGACTTTGTCT
GTTTTACTCTAGTCTTTAGTACTTCAAAATTAACCTTTTTCATATCCATGATCTTGAGTCCATTTGGGGGA
TTTTTAAGAATTTGATGATTTCAATACACTGTTCAAAATTAATTTGTTAATTTTATGATGAGTATG
TATGTTCTGAAGGTGGTCTATTTAAATTTAAACTATTGTAACCTT
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001024024.2
Summary:	This gene encodes a member of the GTP cyclohydrolase family. The encoded protein is the first and rate-limiting enzyme in tetrahydrobiopterin (BH4) biosynthesis, catalyzing the conversion of GTP into 7,8-dihydroneopterin triphosphate. BH4 is an essential cofactor required by aromatic amino acid hydroxylases as well as nitric oxide synthases. Mutations in this gene are associated with malignant hyperphenylalaninemia and dopa-responsive dystonia. Several alternatively spliced transcript variants encoding different isoforms have been described; however, not all variants give rise to a functional enzyme. [provided by RefSeq, Jul 2008]
Locus ID:	2643
MW:	41.6