

Product datasheet for **SC212327**

KCNC4 (NM_001039574) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: KCNC4 (NM_001039574) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: KCNC4
Synonyms: C1orf30; HKSHIIC; KSHIIC; KV3.4
ACCN: NM_001039574
Insert Size: 1102 bp
Insert Sequence: >SC212327 3'UTR clone of NM_001039574
The sequence shown below is from the reference sequence of NM_001039574. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCCCAGGCTGAAGTCTCACCCCTCTCTAAAGCGGCACCAACGTGAGAGAGACAGGCAGACAGACAGAA
AGCCAGAGGCTTAGGGAAACTCTGGAACCCAGACAAGAATCTTTTCGCTGGGAAAGACTCAGATATCCT
TGTTTGACACAGGACTGGTGGAAAATCTCCCATGCGACCCTCGGGGCCAGAGCCATCTGGGTCTGATGT
TCTGTTCCATTGTACATCGAAGAGATATATATGCACATATAGTATCTATATTCATACATATTATACTCT
TGTGTGTAGTGCACGTGCTATTGGTGGTTGTCTTCTTTGTTAGGCTGTGTCTCCCTAAGCCCTTGCCC
CACCCAGAGTTTCCCGTCCCCTTCACTGATTTCTGTTGTTTCTGCTGACTGTGTGGGTGGAATGTCCCA
AGAAAAGTGCATCTGGGAATTGCCAGTCCAGCTGGGTAGTCCCAGGCTCCTGTCTTGGGGATGTTTCCC
CTGTCAGCAAGTAACCTGGTGAAGTCTATTGAAGGCCAGACTGCCCCCTAGGGTCACTGCTTCACTAGC
CGCACCCACCCAGATTGGGGTCTACCTCCCACCCACATCCTCGTTGTGGGGGACTTCCAGGGGC
TCCTCTGCAGCCTCCTCCACTACTTCTCCACCCCATCTATGTCCTTGACTTAGGGGGGACTTTTGTCT
TTTTTAGATTTGATTTTGTCTCTCCTTTGTCTGTTTGTGCAAAGATGCTGCTGGGCAGACAGGC
AGGGAAAGGATCTGTCTGCCATCTGGCCAGGGGTCCGAGAAGGAAGCCTTGGGCAAGAGGAGACC
AGTTGCAATACTGTACTTCTGGTCAAGTGGCCAGAGGATGCGTGAATAGCAGAGGCCAGGTGACCCCT
TCAGCCTTGGCCTCTGCCCTCCCTTGGCCCTCCCTCCCTGCTCCTCCCTGGTGTGGTCACTCCTTTT
CTAAAGCTGTCCCCTCGTGTGTCTGTTGGGATGCCCAGGCTGGGCCCTGTGCCCTGTCTGCATGCCCTC
CAACTGTCATGCTGTGCTCGAGCCCAATAAGACATCTGGAGCATCCTGCTGCTCCTGCTGTGTGA
ACGGGTAAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001039574.3</u>
Summary:	The Shaker gene family of Drosophila encodes components of voltage-gated potassium channels and is comprised of four subfamilies. Based on sequence similarity, this gene is similar to the Shaw subfamily. The protein encoded by this gene belongs to the delayed rectifier class of channel proteins and is an integral membrane protein that mediates the voltage-dependent potassium ion permeability of excitable membranes. It generates atypical voltage-dependent transient current that may be important for neuronal excitability. Multiple transcript variants have been found for this gene. [provided by RefSeq, Jul 2010]
Locus ID:	3749
MW:	39.7