

Product datasheet for SC212328

MTFR1L (NM_001099627) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MTFR1L (NM_001099627) Human 3' UTR Clone
Symbol:	MTFR1L
Synonyms:	FAM54B; HYST1888; MST116; MSTP116
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001099627
Insert Size:	1094 bp
Insert Sequence:	>SC212328 3'UTR clone of NM_001099627 The sequence shown below is from the reference sequence of NM_001099627. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TGTCCAAGGCCAGCAGCTTTGCAGACAATGGGTATCCTGAAGGACTTTCACCGAATGAACAGAGTC
AAGATCTGAACCGGAGTTTATTGAAGGAGGAAGACCCTGCTGTGCTTATCTCTGAGGTCCTAAGGAGGA
AGTTTGTCTAAAGGAAGAAGATATCAGTAGAAAAGGAAATTGACAACCCTCAGCTCTGCAAACTCAGT
CTCATGCTCCTGGAATACCTTCAATAGCTGCCTTCCTCACGCGAGATGTTTCTGCCTCTTAAGGATAGA
TCTTCTGCAACAGTCTTGCTGACAAGCTAGAGCTTGGACTGAAAGAGAAGAGCTGGATTATATATTCC
CAGACTTCAAACCTAGCAGAAGCTAAGGCTTGTGATTTGACCTGAGACATTTGTTTCAGGTAATCGTG
TAGAATGAAGTATCTTAGTTTAAAGGTAAGAGAGAAGTTGTTTCTGGTTTTCTTGCCCTGTGTGA
AAATAGGTCCTAAATGACTGACTTCACTGCATTAGACCTATAGCTGGTCTCACAAGACACTTTGTGCC
CAGCTGTCACTCACTCAGCAGCTTCCTTGACGAGAGCAGGGCTGAGGGGAAGGGGCTATGAATGTTTG
TATACATGTTACAGGGCACGGAAAATCTTATGCTGCTCCGTCATAAACCTACACCAATGCCAGCAAT
CACCTCCTCACTTCCTTGCTAGATGTAGAGGTGAGGCTGCTGAACAGCCAACACATGGGCTACTGC
TGGGAAGCCTGGGCTGTTTTTTTCTTAAACACATTTTATATTACTGAACAACCAATCTACCTCCAC
GGCCTGAGGCCTTATCAGTTCCTAGTATTAATACTTTCTCTTCCACGGACTTTAAGCCCGGTAGGAA
AGAGAGAGGAGGAGGGGAAAGAGCAAACCATCTTCTTCCAGGCCCTTGACTGCTCCTTTGGGCTGGG
CCAAGGTTTGTATGTACCACACCATGCATGACTCAGATGCCCTCAGGTCCTTTCTCTATGGTATGTAT
ACTGCTTGTGTTGGGTTGAAGCACTACCTGACATTAAGGAAGGACTTGAGAGAGAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001099627.2</u>
Locus ID:	56181
MW:	41.8