

Product datasheet for **SC212331**

Phospholipase c eta 1 (PLCH1) (NM_001130960) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Phospholipase c eta 1 (PLCH1) (NM_001130960) Human 3' UTR Clone
Symbol:	Phospholipase c eta 1
Synonyms:	PLC-L3; PLC eta 1; PLCL3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001130960
Insert Size:	1105 bp
Insert Sequence:	>SC212331 3' UTR clone of NM_001130960 The sequence shown below is from the reference sequence of NM_001130960. The complete sequence of this clone may contain minor differences, such as SNPs. Red =Cloning site Blue =Stop Codon

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

CAAGCAGTGATGATAAACAGAAATTTATTTCTTTTGAGACTGTGAATTATTTAAACTGCATTCTTAA
GGTTTTCAAGGATTCTGTAGAATGTCAGAGTTTTAGTAAATAAGTTCCTTGGCTTTGGAATGATTTTC
AAATGTATTTTAACTTGTATTTGGGTCTCCCTTTGACTTCATGTGTTCTCCCTGTTAATGCGTTTGTA
TATAGATTTGGTGACATTTCCCATGTACCTATTATGATCTTTCTCCCTTGAGTTGTTGAAACTGTATCA
GTATGAAATATGTGCATGCCCTAGATGTGAAACTCCTCAGCAGCTCGGGTTATAATGTGTGTGTGTTTCA
CCAAGATGTCCTTTACTGTGTCCTCAGATTGCAGTCACAAATGTAGTCATTTTTACTCTCTAACTGTTGA
GTTGATTAAGAGACTTTGCAAACGACAAAGGGATGAATTCCTTACTTTAATCTGTTATCATTTTTCTAT
GTTTCCCTCTTTGTTGAGAACCCAGATGCATTTTTATAACTCAGTTTTAAAACTTTAAATAGTTACC
TTGCCTTTTAGGATGTTCTTATCCCACCATAATGAGAGTTGAAAGGGGATGGATAGCTGCTCCCATGC
CCTTCCCACTTTTTGGAATAGGCCGTGAGGGTGTGAGGAAGAAGGCTGTCTTTGTACATAAGGACAAAA
TTGTTTGTTTTACATAAATTTTGTACATATTTTGTCTAATGGCTTTGTATGTAACAAGAAGCGAGTTGC
CAAACCTGTTGACTTTTGAATTTTCTGATTGAATTACAGACTGCGAACAACGGCTTTGAGAAATGAG
GGACTTCCATCAGACTCTAATGATAATAGTAGCACAAATTGAAAACCTCCCAAAGCTTTTACAGAAATAT
TTTCTCATAATAAATCCAAGTGAACAGATAATTAGAAGAAACCTTTTCTTCCAGGGAACCAAGCAACT
CTATTTTAGTACTGACATGCATTATTTTCACTGTGAATTCACCTTTTTTATTGCATGTTGAGATGTCCTC
TTTGTTTTTTTTTTTGTAACATTAAGTGAATGATGTTCTTCTGGAATTCATG

ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCG



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001130960.1</u>
Summary:	PLCH1 is a member of the PLC-eta family of the phosphoinositide-specific phospholipase C (PLC) superfamily of enzymes that cleave phosphatidylinositol 4,5-bisphosphate (PtdIns(4,5)P ₂) to generate second messengers inositol 1,4,5-trisphosphate (IP ₃) and diacylglycerol (DAG) (Hwang et al., 2005 [PubMed 15702972]).[supplied by OMIM, Jun 2009]
Locus ID:	23007