

Product datasheet for **SC212340**

FBXO34 (NM_017943) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FBXO34 (NM_017943) Human 3' UTR Clone
Symbol:	FBXO34
Synonyms:	CGI-301; Fbx34
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017943
Insert Size:	1115 bp
Insert Sequence:	>SC212340 3'UTR clone of NM_017943 The sequence shown below is from the reference sequence of NM_017943. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
AAAGGGACTGAAGCTGAAGAGGAATACTAAAGTCCATGTGAGAGGCAACAAAAGGACCGTTTCTAAAG
CTGCAAAACACCTAGATACACCGTTCAAATGAGCGTAGCCCCCTGAGTCATCACTCTAGAAGAATCTGT
ACATCATCAGGACTGCATTGCTCAGGCATTTTCTAACTCTAAATTTACGAGCTGTACAAAAAATTGG
TCTTGTTGTTTATAGTGGCATCTCATGTTGAACCCGGGTGGTATCCACAGTTGGATTCAAGTTGGCTG
TGAATAACTGCCTGTTTTCTAAATCAAACCCATCCTCAAAGGATGAAGACTCACCACCATCCAGGACA
TTCAGAAGAGTTCACTGCAGATGCTGCAGGTAGTCTCAAATGGGTTCCAGAAATGTTTTGAGCACT
GGCAACATATTTGAAATAAGTGAATATTGCTGTGAAAAGAATAGCAGGACTTTTAGATGAAAAGTAT
TCTTAAAAAGAAAAGTCAGGCACCCACCTTAGACCTCGTATGCTTGATCCTGTGAGATTGATGTTTGT
GGCTGGAGGTGGATTTTCATGCCCTGTGGTGTACAGTGTATATAATGTTGTGTTTTCATGGGGCTAT
GAAAGTGCACGTTAACTGAGCGCCTTACCTTTAGATGAGTGCTTTGGCCCTCTGTGAATAGCAGC
ATTAATCCAGTTGTATATAATGGACAGCTAACGGAACAATATAATCACCACAATGCAGCTAGGATAG
TGTTGCGGCTATAATTTGTGTTTTTTTTTTTAAATTGTCTAGTCTTAAATTTGTACATCTTGATATAA
ATATGAATGTTTTCCAAATAAACTATGAATGTTTCTGTATAATATATGAATGTTTCTGAGAAGAACT
CTAAATAGTTGAAAGGCTAACCTGCTCAAAGGATACCAATAATGGTTTAACTGGACAACCTGAAAATT
AGCATAGAAAACAATCCTTTGTTATATTTAGTGATCCACAAGATTGAGAAAATATTATATAGTTAGAT
AATAACATTCTGTCTACTTTATCCTGTCTGGTTACAAAATTTTTTAACTTAAATAAAAAACATGCAT
CTTAAATGGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_017943.4</u>
Summary:	Members of the F-box protein family, such as FBXO34, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]
Locus ID:	55030
MW:	43