

Product datasheet for SC212348

LENG8 (NM_052925) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LENG8 (NM_052925) Human 3' UTR Clone
Symbol:	LENG8
Synonyms:	pp13842
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_052925
Insert Size:	1118 bp
Insert Sequence:	>SC212348 3'UTR clone of NM_052925 The sequence shown below is from the reference sequence of NM_052925. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CTCAGCCTGGCGCAGCTGTCAGCCTTCAGCACCAGCGAGGAGGGGCGGGGCGAGGGGCTGCAGCCC
CCAGCGCTGCCTTTGCGGATTCTGTTTTGAGCCGTGGACTTGGGTTGTAAATTTATTTGTGGGAGTG
CGCTCCAGGAAGAGCCACCATCCCTGCCCGCTTTCCACCGGGGAGTCTGTACAGAGATTTTCTAC
GTTTTTATTTTTGCCTCAGAGGATGGGATTGGGGAGGAGGGATGGGCAGCGAGGGTTGGGGCAT
GGTCTGCAGGCTCATCTGTGTCCGCTTCACTCCAATAATGCTGTCTCAGTGTTTTCTCTCTCTCT
TTCGAGCTTGACTCCGGTACCCGACCCGGCGCCTGGCCCATCCCATGCCGGGGGGCCAGTGGAAGA
AGACAGGCCGTCCAGCCCGTGCCCGCTGCGGCGGGGGACCCAGCAAGCCCGCCACCGCCCGTGCC
TCACCTGCTTCGCCACAGACTCTTGTTCCAGCCCTTGGGGCTCCGTGTTGGGGTGGGGAGCTGC
TTAGAGACTGTGCCCCTCGGCCCCCACCCTGAAGTGCCAGCACCAGCACCAGATCCTCCGCC
GCCACCCGCACTGAGGACAGCCGGCGGGCGCCTCGTCTCAAGTTGTATAAAGTTGTCTCCGTGTC
CCCTCCTCCCTCTGCCCCAGTGTTCTTCTGATTTTTTTCCCTTTCCCTCCCTCCTCTCCGCAT
TCTTCCCTTGGTTCAGCACAGTAAACGGTCCCTCCCTCCCTGCCTTCATGGATCACCAGCTCAG
TCATGTTGCCCTTCTTTTTCTTTGTGTGTGTTATTTAAGTTATTTTCTTCCCTCCTCCTCTTTT
TTTTGGCCCTCCCTCCCTCCTTCTGCCATGTAAGTGGAGGATGTGCTATGAGTTTGAAACAGCT
GGACTGTCAGGCTGCTTTTTTCCAGATGTTCTCCTCTGCTCCCTTCCCTCCTCTCCCTCCTTT
TCCTTCTTCTTCTTCTTCTTGGAGCACTGAGCACCATTGGAAGCTTGAGAGAAACAAAATTAAG
AGAGAAAGAGAGAG
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_052925.4</u>
Locus ID:	114823
MW:	40.1