

Product datasheet for SC212353

WDR92 (NM_138458) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: WDR92 (NM_138458) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: WDR92
ACCN: NM_138458
Insert Size: 1463 bp
Insert Sequence: >SC212353 3'UTR clone of NM_138458

The sequence shown below is from the reference sequence of NM_138458. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
CTGATCGTTACAAAGCTCAATAAAATTTGACTTCAAGACTTTGGACTTGAGGCTGGGCACGGTGGCTCA
CTCCTGTAATCCCAGCACTTTGGGAGGCCGAGGCGGGTGGATCAGAAGGTCAGGAGATCGAGACCATCC
TGGCTAACATGGTGAAACCCCGTCTACTATAAAATATAAAAAAAAAATTTGCTGGGTGTCATGGCCGGC
ATCTGTAGTCTAGTACTTTGGGAGGCTGAGGCAGGAGAATGGCGTGAACCCGGGAGGCGGAGGTTGCA
GTGAGCCGAGATTGTGCCACTGCACTCCAGCCTGGGCAACTGAGCAAGACTGTCTCAAAAAAAAAAAAA
AAAAAAAAAGACTTTGGACTTGAAACCTTCAGAGGAACACTCACACAATTTAAGTAGGTTCTGGAGTC
CTCTGACCCTCATTGAACAAAGCTGGATTTGGCTGATGATAAAAGACCCAGCCACAAGCCTCTGGGCC
AGCCTCCAGCCTCTCTCGGGCTTGTGTGCTGACTTGTGCTGCACGTGGGTATCAGAAGCACCAAGTA
AACTATCTCCAGGGCCAGTTCTTTTCTCTCCCTTCCATTTCTAAGTCAATAGTTGTCTATAATTAACA
TATTTTCAGATTTGTTCTGTGTCCTGAAGTTCCATATTTTATCTTTATGCTTGTGTTATGTTTTAGTGT
CTATAATGTTGTGGTGGGAATCATTTTCAATTTGGTTTAAAAGAGTTTGCCTTGAAAAACTGAATTAACC
TTCATAAATACATTTTGGAGATTTTCAGTGTGAGTAATTTTCTGGACATAATATTATTGTATAAAATAT
CAAAATGGTGGCCAGGCGTGGTAGCTCACGCCTGTAATTCAGCACTTTGGGAGGCTGAGGCGGGTGA
ATCACGAGGTCAGGAATTCGAGACCAGCCTGGCCAACATGGTCAAACCCGTCTCTACTAAAAATACAA
AAAATTAGCTGGGCTTAGTGGCGGGCGCCTGTAATCACAGCTACTTTGGGAGGCTGAGGCAGGAGAATCA
CTTGAACCCAGGAGGTGGAGATTGCAGTGAAGTGGCCACTGCACTCCAGCCCTGGCGACAGA
GTGAGACTCCGTCTCAAAAAAAAAAAAAACAAATATATATATATATCTCTCTCTCAAAACGGCAC
GGTTGGTCTAGATTTGGGTAGCTATTTTCTGCAATTTAATTCATAAAACAAAATGCTTTGATTTTAAAG
AGCTCTTTTTTACTTAAATACAGGAAGTGAATTTGGTGAAGTCTTTGGTTTTTACTGTGGTAAAGTAT
GGTGAATGTACAAGATTGTAGGGTTTTGTAGTCTTATTTTTGGTCATACTTTAAAAATGCATTTTAT
GAATAAATTCTGTTTTTGGATCTACAGAAGCTTTGATTCTGTTACCTAACAGTAGTAATAAAGTTA
ACATTGTCACATCA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```



[View online >](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_138458.4
Summary:	This gene encodes a protein with two WD40 repeat domains thought to be involved in an apoptosis via activation of caspase-3. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]
Locus ID:	116143
MW:	55