

Product datasheet for **SC212370**

ACTR10 (NM_018477) Human 3' UTR Clone

Product data:

| | |
|------------------|---|
| Product Type: | 3' UTR Clones |
| Product Name: | ACTR10 (NM_018477) Human 3' UTR Clone |
| Vector: | pMirTarget (PS100062) |
| Symbol: | ACTR10 |
| Synonyms: | ACTR11; Arp10; Arp11; HARP11 |
| ACCN: | NM_018477 |
| Insert Size: | 1126 bp |
| Insert Sequence: | >SC212370 3' UTR clone of NM_018477 The sequence shown below is from the reference sequence of NM_018477. The complete sequence of this clone may contain minor differences, such as SNPs. Red =Cloning site Blue =Stop Codon |

CAATTGGCAGAGCTCAGAATTCAAGCGATCGC

CTCAACCACCTCTGATGAAGAGAGCATTTCCTACTGAGAAA**TAG**AAGTTTGATTA AAAATCAACCTTGCT
TCATATCAAATATTTAACCAATTATAAGCAAATGTACAAAGTATGTAGGATGTTTTGTTATAGAGGACT
ATAGTGAAGTAAAGCATTCTGTGTTACTCTTGCATTAATATATAATTCTTTTGACTTTGTTTCTCT
TGTGTAGTGGTAAATGGTAGCTGGTGCTTATTGAGATTTGCTGTATTTATATCAATAAAGTATAGTAAA
GCAGTTTGATTTTGAAGTTTGTATGTGGCTTTTTTTTTTTTTTTTTTTTTTTGAGACGGAGTCTCGCTCT
GTCACTTAGGCTGGAGTGCAGTGGCACAATCTACTCATTGCAAGCTCCGCCTCCCGGGTTTACGCCAT
TCTGTCTCAGCCTCCTGAGTAGCTGGGACTATAGGCATACGCCACCCCGCCGGCTAATTTTTGTATAT
TTAGTAGAGACGGGGTTTACCATGTTAGCCAGGATGGTCTCAATTTCTGACCTCGTGATCCACCCTCC
TCGGCCTCCCAAAGTGTGGGATTACAGGCGTGAACCACTGTACCAGCCTATGTTGCTTGCTTTCAAAA
AGAACAGTATTTCTCCACAGTTCACATAACCTGTTCTGGGCTCATAATTGTGTTTCCAGCAGTACTCC
ATTTAAAGCATAAATAGTGTTAACCAATTACCAGTATTTCAAATACTAACATGGGTCCATTTAAAAAAA
TATTTATTTTGGTACTTGGGTGCCATATTATAATAATGTTGTTTTCCAGAAATATATTTGATCTTTGG
ACTGAGAATGTTTGGTAACATAGTCAAGCATTGTTAAGCCAAGTGTGGAAATGTTCACTTTTTAACTTT
ACAGTTTTTTAATGAGCAATTTACATTTCTAAGAAAAAGATACTTCATTTTTATATAAGGTTACAAC
TGCTTTATAAAAATGATTACAAATGTCATTAATCTTTGCGTTTTGTGCTTTATATCTTCTCAGGGCAGAG
GTTCTTGAAATATGACTCAATCATTAACTTTTATTTTATTTTTTTCATTTGTTTTAATGAACCTACCTC
TTGTTG

ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: Sgfl-MluI



[View online »](#)

| | |
|------------------------|--|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | <u>NM_018477.2</u> |
| Locus ID: | 55860 |