

Product datasheet for SC212372

FAM72D (NM_207418) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FAM72D (NM_207418) Human 3' UTR Clone
Symbol:	FAM72D
Synonyms:	GCUD2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_207418
Insert Size:	1129 bp
Insert Sequence:	>SC212372 3'UTR clone of NM_207418 The sequence shown below is from the reference sequence of NM_207418. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AATATCTCAGCAGAGGAGTGATTAGATGAATGGAATTATGATATATATGATATACAACTTTTTTCTA
TTTTAAAAATATTAATGGATCAACTTTAAATTTAGTTGCCAGTGATCTTTTTGGAAAACAAAAA
TGGGGCATTGTTGATTTATTTATTTCTGTCTCTAATTAGTTACCTCAGTTTGATTGAAGCCAGTGGA
GTTGTGCTTTTCTCTACTTCTACTTCTCTCTCCCCACCTTTTTCTGCCAGTGAGGTGATTCTTAA
ATTCAGACGGGAAGATTCTTTCACATATCACTCAGTTACCTCCCAATCTGGGGGAGTTTTCTTACAAC
TTGATACCAAGATACCATTAATTTACATTCCTGAATAAAGGCCTAGTACCCACGCATATTCAACCATG
CATATATCAAGTTCAACCGAGTTTAAATAGGGGATTAAAAAACAAGCTGTTAGGTTTCCATGGGCACT
GGTTCTCATAGTTCTATTGGTGATAACTGCTTTAACATGGAGCAAGAGTTTGTGAATCAGGAAATAGA
ATAAAATTTAAATATATAGAGGAATCCTCTTGATTGCTCAGCATGATGTTAGATAAATGAGTT
TGTCAGAAAATATCAGTATACGCTGTTACCAATGTTATTTATTTACATCTTCTAAAGCCATTATGGA
TATTGTATTATGAGAGCTAAACCTAAATAAGTTATCCTGTTCCCTAGGACCTTCTGTAAATAGTGAA
TTTTAGACGAGTGGTCTGCTCCTAAATCTAAATAGAAAAAACTAAAGCGATTGCTTAAGCCATTG
TACATTATAAAGAGCTGTTTTGTTTGCTTTGCTTTGCTTTGTTTTGTTTTTAAAGCTGCATTGAGA
GCTACAAAGGAATAGGAAAGTAGGGTAGTGTGGATTCTGGTTTTATGTAACCTCTAAATAAATGTATC
TCTTTAATATCTCAGTTGTAGGGATTTGTCAATACCAAGCAGACTGAGTTGTGGTTTTGTAAATAAA
GTTTTTCTAAAAATGACCATTCTTCTTTAATTTTTGTTATGCCACATATTGTATGAAAAATATA
AATAAATAGTACTTAAAGTATAGAA
ACGCGT AAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_207418.3</u>
Locus ID:	728833
MW:	43.7