

Product datasheet for SC212373

FAM72B (NM_001100910) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FAM72B (NM_001100910) Human 3' UTR Clone
Symbol:	FAM72B
Synonyms:	p17
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001100910
Insert Size:	1125 bp
Insert Sequence:	>SC212373 3'UTR clone of NM_001100910 The sequence shown below is from the reference sequence of NM_001100910. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AATATCTCAGCAGAGGAGTGATTAGATAAATGGAATTATGATATATATGATATACAACTTTTTTCTA
TTTTAAAAATATTAATGGATCAACTTTAAATTTGTTAGTTGCCAGTGATCTTTTTTGGAAAACAAAAA
TGGGGCATTGTTGATTTATTTATTTCCGTCTCTAATTAGTTACCTCAGTTTGATTGAAGCCAGTGGA
GTTGTGCTTTTCTCTACTTCTACTTCTCTCCCCACCTTTTTCTGCCAGTGAGGTGATTCTTAA
ATTCAGACGGGAAGATTCTTTCACATATCACTCAGTTACCTCCCAATCTGGGGGAGTTTTTCTTACAAC
TTGATACCCAGATACCATTAATTTACATTCTGAATAAAGGCCTAGTACCCACGCATATTTCAACCATG
CATATATCAAGTTCAACTGAGTTTAAATAGGGGATTAAAAAACAAGCTGTTAGGTTTCCATGGGCACT
GGTTCTCATAGTTCTATTGGTGATAACTGCTTTAACATGGAGCAAGAGTTTGTGAATCAGGAAATAGA
ATAAAATTTAAATTTAAATATATAGAGGAATCCTCTTGATTGCTCAGCATGATGTTAGATAAATGAGTT
TGTCAGAAAATATCAGTATACGCTGTTTACCAATGTTATTTTACATTCTTCTAAAGCCATTATGGA
TATTGTATTATGAGAGCTAAACCTAAATAAGTTATCCTGTTCCCTAGGACCTTCTCTGTAATAGTGAA
TTTTAGACGAGTAGTCTGTCCTAAATCTAAATAGAAAAAACTAAAGCGATTTGCTTAAAGCCATTGT
ACATTATAAAGAGCTGTTTTGTTTTGCTTTGCTTTGTTTTGTTTTGTTTTTAAAGCTGCATTACAGAG
CCACAAAGGAATAGGAAAGTAGGTAGTGTGGATTCTGGTTTTATGTAAGTCTAAATAAATGTATCT
CTTTAATATCTCAGTTGTAGGGATTTGTCAATACCAAGCAGACTGAGTTGTGGTTTTGTAAATAAAG
TTTTTTCTAAAAATGACCATTCTTCTTTAATTTTTTGTATGCCACATATTGTATGTAATAATATAA
ATAAATAGTACTTAAAGTATA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001100910.2</u>
Locus ID:	653820
MW:	43.5