

Product datasheet for **SC212374**

TMEM117 (NM_032256) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM117 (NM_032256) Human 3' UTR Clone
Symbol:	TMEM117
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032256
Insert Size:	1116 bp
Insert Sequence:	>SC212374 3'UTR clone of NM_032256 The sequence shown below is from the reference sequence of NM_032256. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACGACTTCTAAAAGTACACCTACGAACAGACTCGGAGATAGACTTGGAGATAACACAAAAAGCAACCT
TGAGTGTAACTTTAAAAATTTAGTCTTTCTTTGTATATGTAAGGTTTACGTAGTGTTAGGTAAAAAT
ATGAACAATGCCACAACGGTGCTCAACATGCTTTTCTAGGATTCATTGTTTTCTATTTGTATTATAAT
ACACGTGCCTACTGTATACTCAACAGTCCTCTAGAGATTGCTTTTCACAATTGCACAAGCTATTACTGA
CTTTACAGCATAGTGGAAGATTAGCTGATGACCCATGTATCTGATGTTCAACCATAGTGGTGCCTTGAG
ACATTAACTGTTTTTAACTGTACCAGAAATGAAGTGTGGAACAGTTACCTAACCTATTTACATGGGC
GTTTTGTATACTAATTTTGTATCTACACTTGATGTCTGAGCAGAAAACAGAAATAGCTAAATGTGACT
CAGGAAGTATCTCTTGGTTTCTTATTCAGCAGCAGAGTTGGTGACTTTGACAACTGGACTGCAGAGAAA
CATGGTGATCACCTTTTAAATTTTATTGGCTGTCTGCCAAATATAAATACAGATGCAAAATTCAGTAAT
AGGAGATCCATAACCAACATGGGTCACTACTCGTGAAATGTGACTTTCTCCACCAGTAATTGAAATG
AGGTGATGATACCTAATTATGTTTTCTAATTAAGATAAATGCTACTTGATTAATAATCCTGCCCTT
CACCTTTGGGAACAAAGGTTAAGAGACACAGTTGGGCGAACTCTCAAATTTATTGGCATTACACAAAG
TCCCAGACAACCAAGGAAGTGAAGTTTTATCATATGAGAGCAGCACATCCACCATTTACAATATTCG
TATATCTTTCTGCAAAATATGGCTCTGGATAGTGAAATGAAAAACATATGCCAACCTGAGCAAGGGA
ACTCCTCAAAAAATCATGCAGCGGAACCTGTGCAGGTAGAGAAGCCGTGCATGAAAGAAATTTGTTTAAT
GTCTTGTTTTGCGTATGTGTTTTTGTGTTTTTGAAGAACTAAATATTGCACATTAATAAATAAG
AATTATACAGCA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_032256.3
Summary:	Involved in endoplasmic reticulum (ER) stress-induced cell death pathway.[UniProtKB/Swiss-Prot Function]
Locus ID:	84216
MW:	43.7