

Product datasheet for SC212388

Renalase (RNLS) (NM_018363) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Renalase (RNLS) (NM_018363) Human 3' UTR Clone
Symbol:	Renalase
Synonyms:	C10orf59; RENALASE
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018363
Insert Size:	1129 bp
Insert Sequence:	>SC212388 3'UTR clone of NM_018363 The sequence shown below is from the reference sequence of NM_018363. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGGATGATGGCGATTGGATTTCCATCTTGACTTCCTGGAATTGGAGCACACAGTCAGGTTTTATTGTA
TTTTTTTTTTAAGGATACCACTCACAGCCTTAGGATAGCTATTATTTAGAAGCAAAACAGAAGATA
AATGTTGGCAAGGATGTGGAGATATTGGATTCCCTTGTGCAGTGCCGGTGGGAATGTAATGATGTAG
CTACTATGGAAAATGATACGGCAATTTCTTTAGAAATGAAATATAGAATTGCTGTATGATCTGCAGTTC
CACATCTGGATATCTATCCTAAAAGAAAGTAAAGTAGGGACTTGAACGAACATTTGTACCAATGTTCA
CAGCGGCTTTATTCACAACAGCCAAAAGGTGGAAGCAACCCAGTGTCCATGGATAGATGAATAGATAAA
TAAATGTGGTATAAACATACAATGGGCTATTGTTTAGCCTTAAAAGGGAAGGAAATCTGCATGCTG
CAATATGGATGAAGCTTAAAGTCATTATGCAAAGTGGAAATAGCCTATCACAAAAATAATATTACATA
ATTCTACTTATATGAGGAATCTAGAGCAGTCAGTTTCACAGACAGAAAATAGAATGGTGGTTGCCAA
GGGCTGGGAGAAGAGGGCAATGGAGAGTGAGTGTTTAGTGGGTCAGAGTTTTAGTTTGGGAAGGTA
AGTTCTGGAGATGGATGATGGTTATGGGTGCTCAACAGTGTGAATGTACTTAATGCCACAGAAGTGCAC
ATTTAAATGTGGTTAAATCATCACTTTTATGTTATGTATATTTACCACAATAAAATAAGAAAGTTGATA
TTTCTTATACTTACAAAAGAGGAGAAGGGCATTGCAAATCAACAAGAAGTGTGAGGCCCTCTCTAG
CAGAAAAATAGACTAAATCTATTTCTTTATCTTTTAAACATCCTGTTTAAAGGAAAATGCCAAAAAATG
GGAAAAATACACACACAAAATATATATGAACATGTTTTGCCTCATGAGTAATCAAAATGTGTACATA
TGTATGTTTATGTATGTGTTTATATTTAAAATCGTGTCTGCCTTATGAGTAAACAAAAGTATACA
AATTAATAACTATAATGAAACGTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_018363.4
Summary:	Renalase is a flavin adenine dinucleotide-dependent amine oxidase that is secreted into the blood from the kidney (Xu et al., 2005 [PubMed 15841207]).[supplied by OMIM, Mar 2008]
Locus ID:	55328
MW:	43.9