

Product datasheet for SC212394

FIZ1 (NM_032836) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: FIZ1 (NM_032836) Human 3' UTR Clone
Symbol: FIZ1
Synonyms: ZNF798
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_032836
Insert Size: 1104 bp
Insert Sequence: >SC212394 3'UTR clone of NM_032836

The sequence shown below is from the reference sequence of NM_032836. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CACCTGAAGCTGCACCGGGGCATGGACTGACTGCCAGGCTGCGTGCGCCCTGCCCTCCACCCAGCCTCC
TGGACTCGGCCTGGACAGGGGACCTCGGGACTGCGCGTGAGGCCCGGCCCTCCAAATCCAAATCCAG
ACGCAGGCCCTGAAATGAGGGGACCTGACTGGAGAGGTGGGGGCCACCAAAACCCACAAAGGCCCG
GAGCTGGGGGACCACAAACAAACAGGGTCCTTAGCTGGGGCAGGGGAGCCAAATCTAGGGAGAGACTC
CTGAGCCTGAGGTCCCTGGAATGAGTGTGGGTAGCCGTAAAGTCCCAAGACATGGGGACTTTGCAGTGA
GCAATGGGTCTCCACAAGTACCTCTCATCTTGAGAGCCCTAATACTAAAAGATGGGCACCCACCCCCAC
CAAGGGAAGACTGCCCCATTCCCTGAGAGCCATCATTCTTAACGACCTTGATCTGGAGAATGTGGAGGG
AGCATGTCCCTGAATTTTCTAGATCCCTCCAAATGCCACCCACCAGAGTCACTGGTGACCCAGAAAA
TGGATATAGCCGAAATCTGCCTTTCCCTTTTTCATTCCCTGTGCTGAAAGAGGGACAGGGTAGATGC
CCCCTGCCCTCGAATCCCCCTCCCCGACTGTGGAATGGATCGACCCTAACGATCTTCCCCGCCCAAA
CACTAGAAATAGACTGGCCTGAAATCCCCTTGCCAGTAGAATGGACTGATCTATGTGCACACACCCCA
TCACATGGAATGGGCTGGTCTAGGCTGTGGCCTGCCACCTTCTTAGAGTGAATAGGGCGGACACTCCT
CTCTTTTCTGTAGCGTGTGGACCGGTCCACGCAATTCTTATCCTGTGAACATCATCTGAGTGGGAGGT
GGTGGACACTGGGGTTTCTTCCCTCTCTCCGTAGCATCCGTTGGTCTTTCTCTCCATCTCTGTTGGTT
TGTCTGTCTGTCTTCTCCCAATCCCTAGGGGAAGGGGGCATTGGCTAGGGGGTGCCCTGTGAGC
CTCGACCTTGCCCCCTCGTCCCTCTCCCAAGTGTTCAGGACCCCAATAAACCTTGTCTGTGAGTC
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_032836.3</u>
Summary:	This gene encodes zinc finger protein, which interacts with a receptor tyrosine kinase involved in the regulation of hematopoietic and lymphoid cells. This gene product also interacts with a transcription factor that regulates the expression of rod-specific genes in retina. [provided by RefSeq, Jul 2008]
Locus ID:	84922
MW:	40