

Product datasheet for SC212444

DCAKD (NM_001128631) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DCAKD (NM_001128631) Human 3' UTR Clone
Symbol:	DCAKD
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001128631
Insert Size:	1119 bp
Insert Sequence:	>SC212444 3'UTR clone of NM_001128631 The sequence shown below is from the reference sequence of NM_001128631. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CTCACCCACTACCTTCTGCCTTACGCCTAGTGGGGCACTCAAGGCAGGGAGCCCCAGGCCTCCATCTAT
CTCCTTGGAGGCTGAAGCCAGGTAACACATCCTGTTTCCTCTCAGCCCCCTCCCCACACACACACAC
ACACACACACACACACACACACACACACACACACACACACACGACTCTCCGAATCAGGGCTGGGCCTCTGTG
CCAGGCCATCCCGTCTTTGGAGTGTATCCTATTTGAGATCAAGAAATAGCCTGCTGTCCTTCCCACCTT
CCCTGAGCAGCACCTACAGAACAGCAACCGTGATGGTGTGTGAAAGCAACAGTCAGAGACTGCTATACA
TTTCTGTGGGTGGGGTAGGGAATACTGGCTCTGTGGCCCTCCTCAGAGCTGGGAGGCCAAGGGTTCAC
TGGGCTTGAATCCTCTCTAAGCACCTACCCATTGTGGCCTAAACCTCCTCAGTGGACCACGGGAAATC
ATGTTGGACCAAGGGTTTCCATGTTCAAACAGTCAACAAGGCCCTCCCCTAGGTTGAGATATTTTTAT
ACTGAAGGCCAGGCCCCCCAGTATGGCTTGCTGCCAGCAAGCAGCAAGGATTCTCTCGTGGGTAACGC
TGAGCCAGGCAGGACCAGCCTCTTCTGGGAGACCCCTCTCCCACGTAGGGTTTGTGTAGTGCTCCCA
CATCCTGCTTATTGCCTGCCACCCCTGCTTCTCGCCTGGACCTCTTGGTATTCGTGTACACCATCCTT
GCTGTTTCTTGCCTCTGTGCCTTCAACCTGCTGCTCCTGGCCTGGGATACCTTTTATTTCTTTATTT
CTTCTTCTTTTTTTTTTTTTTTTTTTTTTTTAACTACAAGTCTCAGTTACAGCACTGTCTAAAG
CCCCAGGCTGGGTAGGTGGTCTAGGATACCAACTCCCCCTCAACATTTTACCTTAATCCAACAGTACA
GTAATTAGCCGAACTTGCTGTGTTTTTCTACTAGACTGTAAGCCCTCTAGGGACAGGGACAGTGTCT
TATTCTCTCTGTGTTCTTAGTGAGAGCCTGGCCCGCTGTTGGTGTAAATAAAGTGTGTTGATCTGA
AATGCCCCAGACCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001128631.3</u>
Locus ID:	79877
MW:	40.9