

Product datasheet for SC212467

NIPSNAP2 (NM_001483) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NIPSNAP2 (NM_001483) Human 3' UTR Clone
Symbol:	NIPSNAP2
Synonyms:	GBAS
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001483
Insert Size:	1137 bp
Insert Sequence:	>SC212467 3'UTR clone of NM_001483 The sequence shown below is from the reference sequence of NM_001483. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATCCCACTGAAGACCTCGCCCTCCAGTAAGCTGTAGAGTTTCTATGTGCCTACATACATTTCTGTGA
CAAGTATTTGTCGTAAATTAATTTTAATTGTGTATCAAGTGAAAAAGAACTGAGGTTTTAAGCTGC
TGTATATAGCTTGTGAGAAACCTCTTTTCTTTAAATTTACATAATCACAAGAAAGGAAAGAATTACAG
TTGGACTGATTGTGACAGTGCCTTGTGTCCTCTTTGAAACACCCCGTGTGTCCAGTATACCTTATAA
CACTGAGCCACTTCTCCCAACCTCCAGAAGGGGTCCACGTTGAATTCTGAATCATCTTGAAAAATAAGA
TTCCAACCACAAAAAAATTTAGCCATTTCTTTACTAAAAAAACCAAAAAACAATCTGTTTTATAAT
CACAGATTTTTAGACAAATTTCTTGATCAGGAAGAAATACAAATTTTGTGATGTTTCTCAAGCAGTTT
TTCTGAGTAGTTTCTGAGGAGGAACAAATTACAAGTGACCAATAACTGAAAAATGTTTTAACTCACTC
TCATTTGTAAGCAGTCCACATAGTAGACAATGGGTTTTCCAAGCTGGGCAAGGTACATTTAATCAGTAA
ATCAGTTTCACATCATGTATTGTGATGTTTCAATGTGAGACACAAAAACAATGGCTTGAACTTGTGTA
TCATATGTGATTTTGAATGAACACCTTGAATAGCACTAATTTTTATTTGTGGTATTTTCTATAACAA
AACAAAGTAGCTTAGGAAAAGAGGTTTTATTTGTAAACGATCATTGTGACCTCAGACACTCTCTGGC
TAATATTTTAATAAGCTCACAGCAGATAATTCTGAGATCATGGGTGAGGGGTGGTGCATGTTGAGATTT
AAATTGGCATAAAGCTGCATACTTTTGTCTAGCTGTTTGATTTCATTTTTAATATAGTATGCCAATT
TTGTGACTGTTACCATGTGAAAGTCCTGTTGAAATGAACAATTGTCTGCCCCACAATCAAGAATGTATG
TGTAAGTGTGAATAAATCTCATATCAATGTCAAACCTTTACATGTGAATGATTTTCTCAAAGAACAT
AGAAAAGTCAATAAAATCCTCTTAATTTCCACA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001483.3</u>
Summary:	This gene encodes a member of the NipSnap family of proteins that may be involved in vesicular transport. The encoded protein is localized to mitochondria and plays a role in oxidative phosphorylation. A pseudogene of this gene is located on the long arm of chromosome 2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Feb 2011]
Locus ID:	2631
MW:	44.3