

Product datasheet for **SC212476**

CTNNB1 (NM_001904) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CTNNB1 (NM_001904) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CTNNB1
Synonyms: armadillo; CTNNB; EVR7; MRD19; NEDSDV
ACCN: NM_001904
Insert Size: 1131 bp
Insert Sequence: >SC212476 3'UTR clone of NM_001904
The sequence shown below is from the reference sequence of NM_001904. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGCTGGCCTGGTTTGATACTGACCTGTAAATCATCCTTTAGGTAAGAAGTTTTAAAAAGCCAGTTTGG
GTAAAATACTTTTACTCTGCCTACAGAACTTCAGAAAAGACTTGGTTGGTAGGGTGGGAGTGGTTTAGGC
TATTTGTAAATCTGCCACAAAAACAGGTATATACTTTGAAAGGAGATGTCTTGGAACATTGGAATGTTC
TCAGATTTCTGGTTGTTATGTGATCATGTGTGGAAGTTATTAACCTTAATGTTTTTGGCCACAGCTTTT
GCAACTTAATACTCAAATGAGTAACATTTGCTGTTTTAAACATTAATAGCAGCCTTTCTCTCTTTATAC
AGCTGTATTGTCTGAACTTGCAATTGTGATTGGCCTGTAGAGTTGCTGAGAGGGCTCGAGGGTGGGCTG
GTATCTCAGAAAGTGCCTGACACACTAACCAAGCTGAGTTTCTATGGGAACAATTGAAGTAAACTTTT
TGTCTGGTCCTTTTTGGTCGAGGAGTAACAATACAAATGGATTTGGGAGTGACTCAAGAAGTGAAGA
ATGCACAAGAATGGATCACAAGATGGAATTTATCAAACCCTAGCCTTGCTTGTTAAATTTTTTTTTTTT
TTTTTTAAGAATATCTGTAATGGTACTGACTTTGCTTGCTTTGAAGTAGCTCTTTTTTTTTTTTTTTTT
TTTTTTTTGCAGTAACTGTTTTTAAGTCTCTCGTAGTGTAAAGTTATAGTGAATACTGCTACAGCAA
TTTCTAATTTTTAAGAATTGAGTAATGGTGTAGAACAATAATCATAATCACTCTAATTAATTGTATC
TGAATAAAGTGTAACAATTGTGTAGCCTTTTTGTATAAAATAGACAAATAGAAAATGGTCCAATTAGTT
TCCTTTTTAATATGCTTAAAATAAGCAGGTGGATCTATTTTCATGTTTTTGATCAAAAATTTGGGAT
ATGTATGGGTAGGGTAAATCAGTAAGAGGTGTTATTTGGAACCTGTTTTGGACAGTTTACCAGTTGCC
TTTTATCCCAAAGTTGTTGTAACCTGCTGTGATACGATGCTTCAAGAGAAAATGCGGTTATAAAAAATG
GTTCAGAATTAACCTTTAATTCATTC
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001904.4</u>
Summary:	The protein encoded by this gene is part of a complex of proteins that constitute adherens junctions (AJs). AJs are necessary for the creation and maintenance of epithelial cell layers by regulating cell growth and adhesion between cells. The encoded protein also anchors the actin cytoskeleton and may be responsible for transmitting the contact inhibition signal that causes cells to stop dividing once the epithelial sheet is complete. Finally, this protein binds to the product of the APC gene, which is mutated in adenomatous polyposis of the colon. Mutations in this gene are a cause of colorectal cancer (CRC), pilomatrixoma (PTR), medulloblastoma (MDB), and ovarian cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2016]
Locus ID:	1499
MW:	43.9