

Product datasheet for **SC212484**

RBM11 (NM_144770) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RBM11 (NM_144770) Human 3' UTR Clone
Symbol:	RBM11
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144770
Insert Size:	1121 bp
Insert Sequence:	>SC212484 3' UTR clone of NM_144770

The sequence shown below is from the reference sequence of NM_144770. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

GTAGCCAAAAGTTCCGAAAGTCTAAGAAGAAGAAAAGATAC**TAG**TATTACCTACAAATGAACTTACCTA
 CACTGATCTTAGTTCTCTTATGAAAAAATAAGATGTTATCCCATCAAATAACAATGTCATGGCTATCT
 TGTCTTTTGAAATATGTAGTAATAAGACTTATCTCTGGTTTGCTGTTCAATTTTGCCTTGAACGACAAA
 AGCTTTTCTAAAAATAGTATAATGTACCACTTTTGTATTGTGCATGATATTTTGGCCATTGGCAACAA
 ATAGACTATTGTCATTTTATTAGTACCAATGATGAAATTTACCAAACTATTAATCAAAAGTCATTATT
 GAACACATTGATAAAGCATCATGTTCCCTTTGATTATCCCTATAACCAACACTCAACTACAAATAGGCT
 GGGATTCAATGTGGAGTGGTGGATTTCTGAAAGAAAGTAGTGCAGATATTTAAGGATTGCATATTGTTT
 CTAATGAGGATGTAATAAACAGGCATTAGTACCCCTTTCTTTAGTACCTAAATATGTTGTTTCTGAGCC
 ATATTCTTCTCAGGTTTATTTTATAAAGTGAGGACATATGATAATGATACCAAGGAAAGATAATTATC
 ATGAAACCAACTTTTGGATCATTCCCTTTTATACAATTGTGTATTTATAGGGATACCTTCTCTCTCA
 CATATTAATTGTATGTCATTAATAATCAGTTGGCTATGTCAGTGTGTATTATTCACCTCTTTGCTAGCAA
 GGTGTGAGATGCTATATTAACATTTAAATCATGAAATGAGCAAGTAAATAGAATAATATTTTATTAATG
 TATGATATGTAATACTTATTGTCTAAATCAAGTTTCCTATTTTGGTTTATTTTATTTAGAGAAAAATA
 AAATTATTTGGGCATTTTGTGTTATAGCTCATTACAACTAGAGTGATGAATATATATTTACTAAAGC
 TAATTGAAAGAGTCATGCAACTCAGTTGGGAATTAACCACTCAATTCAATAGTAAAGCCATGCTTGGAT
 GTTCTTTTGCATGTTGCCTAGCATTTTATAGATTTTATTTAAAAAATAAAGATATGTTTTTAAT
 G

ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCG

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_144770.3
Summary:	Tissue-specific splicing factor with potential implication in the regulation of alternative splicing during neuron and germ cell differentiation. Antagonizes SRSF1-mediated BCL-X splicing. May affect the choice of alternative 5' splice sites by binding to specific sequences in exons and antagonizing the SR protein SRSF1.[UniProtKB/Swiss-Prot Function]
Locus ID:	54033