

Product datasheet for SC212503

ARHGAP36 (NM_144967) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ARHGAP36 (NM_144967) Human 3' UTR Clone
Symbol:	ARHGAP36
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144967
Insert Size:	1130 bp
Insert Sequence:	>SC212503 3'UTR clone of NM_144967 The sequence shown below is from the reference sequence of NM_144967. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
AAACTGGCGTCAGCTACTTCTTCTTCTATAGGTGCCAGACAGGGGAAAAGG
GTGGGGGTACATCTGGGATGTCACAGGAAACATTAAGGAGAGAGTTGAAGGTAAGATCTGAAGGTAAG
AAGGAGTTCACCTGATGCTCGGGTCAGGATGAGAATCCAAACACACTGCCAGCCCTTCACTGGGGA
TGCTTGGTCTCTTCTGCTGGTAAAAGCAGAGATGTTTCTGTGCATGCCAAGCTCCCCGGTGCTACCT
TGCCCTTCTCTTTTACCCTGATCTTGGCTTTCTCTCTCTCTGCAGACTTTCCTTTAATTGATGTGA
CATTTGTGGTAAACACCTTTCCAGGGAACCTCACAATCTTGAGATGCTTCCCTTCCCAGATGGGA
TTGCATGATTTCCCTGACTTTCTACCCTCCTCCAGAGAGCTCAGTTGGAAAGGCCCTCAAGAGGCATG
CTAGAACGTTAGGTCAGCCTACTGACAGCTGACAAACAATTAATGCGAAATCATGTCACACCAACCCAT
AGCCGTGTCCACGCAGCAACTCCACCACCTTAGGATTTCCCCCTCCAAATTATTCAGACCAATGGCTTG
CCAAATGGCCTCTCCAAAATTCTGTACAGTTTGTCTCAGGTACGCCAACAGGGAAACCTCAAGTGTA
GGTCTAATTAGTGTCTTGGGATCCAAAGTTAGAGGAAAATTTAGATTTTATTCCTGGATCTGCTTTA
AAGACAATTGGTGTTTACACCCTCTGTGACGAAAACAGCTAGTTAGGTAAGGACATATAGTTCCAAGT
AGGTAAAGTCACTTGATTACAAATGTTCTTAACATATCGTCTCTGTAATTCCTTTATACAGGACAGTACA
AAATTGTGGGACATGCTCTGGTAACACACAGATATGGGTTGCATATGATCCAGAATTACAGCTGATATT
ATGGATGACAACCTGCTAAGGTCCATAAAATGAAGACTGTATTGTATTGAGGGATAGAAAATTGATCATTT
AATGGGTAACAACCTGCTGAGCTCAAAGATTTGTGATTGTTAAACTTCTCTGGCATTTAATCATTAAATA
AACATCTGTATTGTGACAGCAGCATA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_144967.4</u>
Summary:	GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.[UniProtKB/Swiss-Prot Function]
Locus ID:	158763
MW:	41.7