

Product datasheet for **SC212504**

THSD7B (NM_001080427) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	THSD7B (NM_001080427) Human 3' UTR Clone
Symbol:	THSD7B
Synonyms:	KIAA1679
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001080427
Insert Size:	1142 bp
Insert Sequence:	>SC212504 3'UTR clone of NM_001080427 The sequence shown below is from the reference sequence of NM_001080427. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TTAGCCTACGATGGAGACTTAGACATGTACTGAAAAAGAAATCCAAATGTAGACATCAACTGCCTTA
ACCGCTTTCTCTTTGTAGCTCTCAGACTTCTCAGTTTTTTGAGGAATCTCAAGATGTGATATATTGGG
CAGAATACAAATATTGCAAAAGTAATATTGCCTCAACTTCATTTGGACATGGAGTCAAGGATTATTAGG
TCTGCCATTTTGTTCAGTTGTTGTGGGTGTGTTTTATTTTTTTGGTTTTCCCAAGGGACCTGAAA
ACCCCTTCTCTTCCTGTTTGGAAGTGGGAGGAAGAAAACATGATGGAATCCACAGACTTGAGTAACT
TGATCTTCAGCAGCATAATGACAATCCAGAGGAAAATACAGTCAAGGCATTTACTCTAAATGACGAGTC
TGACACTGCATTGTTACGCACTATATTAGTGCAAAGTCTTTATTCTTCCCATACTTCAACACTGAGTTT
TCTAGAGTTTACTTTGGTTTAAAGACTTTCAAATTGGATTGCCTATTTTTCATGAACACAGAGAGAATG
GATTACCATTTTCAGAAATTCTCTGAGTTTTTAACCTTTAAATATTGTATTTTGTGTTGTAGCCAGGGGA
TGATGGCGCTTCATGGGTTGCAGCTACTGAAAATAGCAGCGTGTGTGTAATTGCTGGACTAGATGAAAG
CTAGGTCATTTCTGAAGGAATGTGTACTGAATGTTAGAGTGACAAATGAAATATGTGGTTAAATTGG
AGAATGAGGTAGATTATTTGATTACTAAACTGTATTTTAACAAAACTTATCCATGTAGATATAGCAT
TAACCACACAGTTGTAAATTCAGTTTAAATGATGACAAACTCTGCTTTTGTAAATTTCAATTTTCTTATC
TGAATATTTATAAATCTTCTTTCAAATTTAATTATCTGACCTCATTTAATATACATCAAACACCGATC
CTGTTTGTACAAAGTCTTGCTTTTATAAGGTTTCAATAATATCTAAAACAACACATTAAGGCTGAGA
CCATTTTATGAAGATAATTGTTTGTAAATCATAGGTGTTGAAAGTAAAAAGGTGCCATCTTGTGGTATTG
ACTTGATTTTATAACAAATAAACTGCTCAAGAGACTGC
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001080427.1</u>
Locus ID:	80731
MW:	45.2