

Product datasheet for **SC212510**

PRC1 (NM_199413) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PRC1 (NM_199413) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PRC1
Synonyms: ASE1
ACCN: NM_199413
Insert Size: 1121 bp
Insert Sequence: >SC212510 3'UTR clone of NM_199413
The sequence shown below is from the reference sequence of NM_199413. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATCCTCAATTCAACCAACATCCAGTCCTGAGAAGCCCTGATCAGTCAACCAGCTGTGGCTTCTGTGCC
TAGACTGGACCTAATTATATGGGGGTGACTTTAGTTTTTCTTCAGCTTAGGCGTGCTTGAACCTTGCC
CAGGTTCCATGACCATGGGCCTAACTTAAAGATGTGAATGAGTGTACAGTTGAAAGCCCATCATAGGT
TTAGTGGTCTAGGAGACTTGGTTTTGACTTATATACATGAAAAGTTTATGGCAAGAAGTCAAATTTT
AGCATATGGGCCTGACTTCTCTACCACATAATTCTACTTGCTGAAGCATGATCAAAGCTTGTTTTATT
TCACCACTGTAGGAAAATGATTGACTATGCCCATCCCTGGGGTAATTTTGGCATGTATACCTGTAAC
AGTAATTAACATCTTTTTTGGTTAGGCATGTTCAATTAATGCTGTAGCTATCATAGCTTTGCTCTTACC
TGAAGCCTTGCCCCACCACACAGGACAGCCTTCCTCCTGAAGAGAATGCTTTGTGTGCCGAAGTTG
AGATGGCCTGCCCTACTGCCAAAGAGGTGACAGGAAGGCTGGGAGCAGCTTTGTTAAATTGTGTTCA
TCTGTTACACAGTGCATTGCCCTTGTGGGGTATGCATGTATGAACACACATGCTTGTGCGAACGCT
TTCTCGGCGTTTGTCCCTTGGCTCTCATCTCCCCATTCTGTGCTACTTTGCCTGAGTCTTCTACC
CCCGCAGTTGCCAGCCACATTGGGAGTCTGTTTGTCCAATGGGTTGAGCTGTCTTTGTCGTGGAGATC
TGGAACTTTGCACATGTCACTACTGGGGAGGTGTTCTGTCTAGCTTCCACGATGAGGCGCCCTTTT
ACCTATCCTCTCAATCACTACTCTTCTTGAAGCACTATTATTTATTCTTCCGCTGTCTGCCTGCAGCAG
TACTACTGTCAACATAGTGTAAATGGTTCTCAAAAGCTTACCAGTGTGGACTTGGTGTTAGCCACGCTG
TTTACTCATACAGTACGTGCTGTTTTTAAAAATACAATTATCTTAAAAATAAATTTAAATCTGTA
TACTTACATTTCAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

| | |
|------------------------|--|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq: | <u>NM_199413.3</u> |
| Summary: | This gene encodes a protein that is involved in cytokinesis. The protein is present at high levels during the S and G2/M phases of mitosis but its levels drop dramatically when the cell exits mitosis and enters the G1 phase. It is located in the nucleus during interphase, becomes associated with mitotic spindles in a highly dynamic manner during mitosis, and localizes to the cell mid-body during cytokinesis. This protein has been shown to be a substrate of several cyclin-dependent kinases (CDKs). It is necessary for polarizing parallel microtubules and concentrating the factors responsible for contractile ring assembly. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2012] |
| Locus ID: | 9055 |
| MW: | 41.4 |